

always by your side




**MEDICAL
INSTRUMENTS**
always by your side



NINGBO SEASIGHT MEDICAL INSTRUMENTS CO., LTD.
Office Add: 9F, Ningshing Building, No.666 Jiangnan Road, Ningbo 315040, China.
Factory Add: No.508, Yindong Road(N), Yinzhou Economic
Development Zone, Ningbo315145, China
Tel: 0086-574-81877115, 89078322 Fax: 0086-574-81877057
<http://www.seasight-med.com>



NINGBO SEASIGHT MEDICAL INSTRUMENTS CO., LTD.



Contents

Breathcare Anesthesia.....	001~016
Urinary.....	017~030
Hypodermic.....	031~047
Examination.....	048~053
Wound Dressing.....	054~067
Medical Non-woven & Accessories.....	068~075
Others	076~084

ABOUT US

OUR COMPANY



Establishment

SEASIGHT is established in the year 2003 as a professional manufacturer and supplier for medical and lab products in China. Our 1st factory- MFLAB Medical with production area 15,500 m2. It includes 100000 grade purification workshops and cleaning rooms, totally 9,200 m2. Our 2nd Factory-SEASIGHT Medical (production area covering 52,000 m2) will open soon.

Team Work

SEASIGHT has a high efficient and professional team, full of wisdom and passion, with more than 20 years of working experience. We take the diligent, professional and responsible spirit, committed to provide customers good products and service.

Export Reputation

Our SEASIGHT products have been exported to over 150 countries worldwide through a network of distributors and subsidiaries, covering North American, South American, European, Asian and African countries.

Leading Products

SEASIGHT has a wide range of products in business with leading products respiratory products, disposable medical products, medical equipment and laboratory products. We can supply clients with OEM or ODM or OBM.

Brands

There are Seven brands in our group: Seasight, Sunflos, Medland Care, Fyocare, Foyomed, FoyoLab and Mflab in distribution.

Quality System

Our factory MFLAB Medical audited by TUV Rheinland every year, more than 40 products has got CE certificates issued by TUV Rheinland, and also passed by GMP, ISO13485, US FDA on-site Audit, Saudi Arabian FDA on-site Audit, Malaysian MOH on-site Audit, and Korea MOH on-site Audit. We have also got a number of patents.

Principle

We SEASIGHT adheres to "Always by your side" principle to strive to achieve excellent quality, competitive price and first-rate service in operation, and cherish every business opportunity.

BREATHCARE ANESTHESIA



LB1110 Simple Oxygen Mask

Made of clear, medical grade PVC
Adjustable nose clip
Anti-crush tube
2m oxygen tubing
The length can be customized

Cat. NO.	Size
LB111000	XS, Neonate
LB111001	S, Pediatric Standard
LB111002	M, Pediatric Elongated
LB111003	L, Adult Asia
LB111004	L+, Adult Europe
LB111005	XL, Adult Elongated



LB1210 Venturi Mask

Made of clear, medical grade PVC
Two color coded diluters for seven concentrations
Green diluter for 24%,26%,28%,30%
White diluter for 35%,40%,50%
Adjustable nose clip
Anti-crush tube

Cat. NO.	Size
LB121000	XS, Neonate
LB121001	S, Pediatric Standard
LB121002	M, Pediatric Elongated
LB121003	L, Adult Asia
LB121004	L+, Adult Europe
LB121005	XL, Adult Elongated



LB1220 Adjustable Venturi Mask

Made of clear, medical grade PVC
Single diluter for one concentration,
seven diluters can be selected
Adjustable nose clip
Adaptor for high humidity
Anti-crush tube
With seven diluters 24% (blue) 28% (yellow)
31% (white) 35% (green) 40% (pink)
50% (orange) 60% (red)

Cat. NO.		Size
Without 60% (red)	With 60% (red)	
LB122000A	LB122000B	XS, Neonate
LB122001A	LB122001B	S, Pediatric Standard
LB122002A	LB122002B	M, Pediatric Elongated
LB122003A	LB122003B	L, Adult Asia
LB122004A	LB122004B	L+, Adult Europe
LB122005A	LB122005B	XL, Adult Elongated



LB1310 Oxygen Mask with Reservoir Bag

Made of clear, medical grade PVC
Standard connector, non-rebreathing
With 1000ml/600ml reservoir bag
Anti-crush tube
2m oxygen tubing
The length can be customized

Cat. NO.	Size
LB131000	XS, Neonate
LB131001	S, Pediatric Standard
LB131002	M, Pediatric Elongated
LB131003	L, Adult Asia
LB131004	L+, Adult Europe
LB131005	XL, Adult Elongated




LB1610
Oxygen Connecting Tube

Made of medical grade PVC
The length of the tube can be customized.



LB1510
Tracheostomy Mask

Specialized mask for tracheostomy
Swivel 360° adaptor with 18F/22M port
Size: M / L



LB141/LB1412/LB1413
Nasal Oxygen Cannula

Made of clear, medical grade PVC
Dipped prong, straight prong, curved prong is available
Color: green transparent, white transparent and light blue transparent is available.
Length: with 2m tubing, different length can be customized, anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Size	Cat.No. : LB1411 (With Curved Prong)	Cat.No.: LB1412 (With Straight Prong)	Cat.No. : LB1413 (With Dipped Prong)
XXS(Neonate)	/	LB141201	LB141301
XS(Infant)	/	LB141202	LB141302
S(Child)	LB141103	LB141203	LB141303
L(Adult)	LB141104	LB141204	LB141304



LB2011
Nebulizer Mask

Operates even at horizontal position Capacity: 6cc,10cc
2m oxygen tubing
Different type of nebulizer chamber are available.

Cat. NO.	Size
LB201100	XS, Neonate
LB201101	S, Pediatric Standard
LB201102	M, Pediatric Elongated
LB201103	L, Adult Asian
LB201104	L+, Adult Europe
LB201105	XL, Adult Elongated



LB2012
Nebulizer Kit

With 1pc adult mask and 1pc child mask
T connector, mouth piece, 2m oxygen tubing
Different types of nebulizer chamber are available.



LB2013
Nebulizer with Mouth Piece

Hand-held, with mouth piece
T connector, mouthpiece, 2m oxygen tubing
Transparent
With/without corrugated tube
Different types of nebulizer chamber are available.

Nebulizer chamber



LB20A

Volume: 6cc
Under the condition of gas flow:
7L/min, atomizing volume is a
200mg/min, in every 200mg,
The percentage of particulate
with the diameter less than 5um
is more than 70%.

LB20B

Volume: 6cc
Under the condition of gas flow:
7L/min, atomizing volume is a
280mg/min, in every 280mg,
The percentage of particulate
with the diameter less than 5um
is more than 70%.

LB20B Patent Number :
ZL 2014 2 0766313.3

LB20C

Volume: 10cc
Under the condition of gas flow:
7L/min, atomizing volume is a
250mg/min, in every 250mg,
The percentage of particulate
with the diameter less than 5um
is more than 70%.

LB20C Patent Number :
ZL 201520637654.5

LB20D

Volume: 6cc
Under the condition of gas flow:
7L/min, atomizing volume is a
about 150-450mg/min,
The percentage of particulate
with the diameter less than 5um
is more than 70%.

LB20D Patent Number :
ZL 201520289256.9



LB20E



LB20F

LB20E

Volume:10cc
Under the condition of gas flow: 10L/min, atomizing volume is 120-450mg/min,
The percentage of particulate with the diameter less than
5um is more than 71%.

LB20F

Volume:10cc
Under the condition of gas flow: 10L/min, atomizing volume is 120-450mg/min,
The percentage of particulate with the diameter less than
5um is more than 71%.



LB60
Suction Connecting Tube

- 1) Made of transparent material for better visualization
- 2) Striated walls of tube provide superior strength and anti-kinking
- 3) Multiple choices of length
- 4) Tube length can be customized

Cat. NO.	Length (cm)	Size (inch)	Type
LB60100118	180	1/4	Long connector
LB60100218	180	1/4	Short connector
LB60100318	180	1/4	Trumpet connector
LB60200118	180	3/16	Long connector
LB60200218	180	3/16	Short connector
LB60200318	180	3/16	Trumpet connector

LB61
Suction Connecting Tube with Yankauer Handle

- 1) Made of transparent material
- 2) Can be combined with different types of yankauer handles and connecting tubes to meet various requirements
- 3) Multiple choices of length
- 4) Can choose different connectors

Cat. NO.	Length (cm)	Size (inch)	Type
LB61110118	180	1/4	Crown tip with vent
LB61120118	180	1/4	Crown tip without vent
LB61130118	180	1/4	Flat tip with vent
LB61140118	180	1/4	Flat tip without vent
LB61210118	180	3/16	Crown tip with vent
LB61220118	180	3/16	Crown tip without vent
LB61230118	180	3/16	Flat tip with vent
LB61240118	180	3/16	Flat tip without vent



LB3011
Anesthesia Mask

Made from Latex/DEHP free PVC material.
Air cushion assures comfortable face fitting.
Tail Check valve/Upright Check valve
Without check valve type are available.
Single Patient use.
Color-Coded hook options.

Cat.No.	Size	Color	Connector
LB301100	size0# - Neonate	White	15mm OD
LB301101	size1# - Infant	Pink	15mm OD
LB301102	size2# - Pediatric	Yellow	22mm ID
LB301103	size 3# - Adult-small	Green	22mm ID
LB301104	size 4# - Adult-medium	Red	22mm ID
LB301105	size 5# - Adult-large	Blue	22mm ID
LB301106	size 6# - Adult-Extra large	Orange	22mm ID



LB3021
Soft Anesthesia Mask

Made from Latex/DEHP free PVC material.
Soft adjustable face cushion assures comfortable face fitting.
With Upright Check valve.
Single Patient use.
No scents or with Vanilla scents are available.

Cat.No.	Size	Color	Connector
LB302100	0# - Neonate	White	15mm OD
LB302101	1# - Pediatric-Small	Pink	15mm OD
LB302102	2# - Pediatric-Medium	Yellow	15mm OD
LB302103	3# - Pediatric-Large	Green	22mm ID
LB302104	4# - Adult-Small	Red	22mm ID
LB302105	5# - Adult-Medium	Blue	22mm ID
LB302106	6# - Adult-Large	Orange	22mm ID



LB3012
PVC free Anesthesia Mask

Made from TPE and PP Material, contain no PVC and are DEHP free
The products provide a comfortable superior airtight fit with a clear shell to enhance patient face visibility.
This products ancillary use with anesthesia breathing circuit.
Non-reusable, single use only

Cat.No.	Size	Color	Connector
LB301201	Size1# - Neonate	Blue	15mm OD
LB301202	Size2# - Infant	Orange	15mm OD
LB301203	Size3# - Pediatric	Pink	22mm ID
LB301204	Size4# - Adult-S	White	22mm ID
LB301205	Size5# - Adult-M	Yellow	22mm ID
LB301206	Size6# - Adult-L	Green	22mm ID



LB3030
Silicone Anesthesia Mask (One-piece)

All-in-one piece soft silicone cushion for better seal and comfort.
Ergonomic shape for lower dead space and better air circulation.
Made of 100% clear, medical grade silicon, latex free
Reusable, all parts are fully autoclavable up to 134°C
Packing: Bulk or PE pack available

Cat.No.	Size	Connector
LB303000	Size0# - Neonate	15mm OD
LB303001	Size1# - Infant	15mm OD
LB303002	Size2# - Pediatric	22mm ID
LB303003	Size3# - Adult-S	22mm ID
LB303004	size4# - Adult-M	22mm ID
LB303005	size5# - Adult-L	22mm ID
LB303006	size6# - Adult-XL	22mm ID

Anesthesia Breathing Circuit

Anesthesia breathing circuit is bridge of anesthesia machine with Patient, input anesthesia gases and output expiratory air from the patient, to accomplish the normal exchange of oxygen and carbon dioxide gas.

Disposable anesthesia circuit can prevent patients from infecting, establish manual airway and match with anesthesia machine and respiratory. Non-toxic, insipid, made of medical degree material, flexible transparent, air tightness.

LB431 Anesthesia Breathing Circuit-normal(Adult/Pediatric)	LB431E Anesthesia Breathing Circuit-extendable(Adult/Pediatric)
LB431R Anesthesia Circuit-Reinforced(Adult/Pediatric)	LB4330 Anesthesia Circuit-Coaxial(Adult)
LB432 Breathing Circuit-Normal(Adult/Pediatric)	LB432R Breathing Circuit-Reinforced(Adult/Pediatric)

LB4410E LB4410 LB4410R double sealed OD: 15mm	LB4610 LB4610S
Catheter Mount (Normal/Extendable/Reinforced)	LB4610 Water Trap

Anesthesia Breathing Bag

0.5L 1L 2L 3L	0.5L 1L 2L 3L
LB4910 Anesthesia Breathing Bag-Latex Free Specification: 0.5L , 1L, 2L, 3L	LB4920 Anesthesia Breathing Bag-Latex Specification: 0.5L , 1L, 2L, 3L

<p>Filter Disposable breathing filter made of medical degree material, mainly used for filtering bacteria and increase gas. humidity degree, effectively prevent patient infecting and mainly used as medical accessories. Efficiency filtration, humidity exchange function.</p>			
LB411001 BV Filter	LB421003 HME Filter	LB421002 HME Filter	LB421001 HME Filter
LB411002 BV Filter	LB421006 HME Filter	LB421004 HME Filter	LB421005 HME Filter



LB3110
Silicone Resuscitator

Built-in 40-60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, silicone mask and oxygen tubing.
100% medical level silicone material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized labels are available.
Built-in 30mm exhale port for PEEP or filter attachment.
All parts (excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C.



LB3110L
Liquid Silicone Resuscitator

Built-in 40-60cm H₂O pressure relief valve
The set includes oxygen reservoir bag, silicone mask and oxygen tubing
100% Medical level silicone material
Latex-free components
All-in-one intake valve connects to oxygen reservoir bag directly
Additional accessories and customized label are available
Built-in 30mm exhale port for PEEP or filter attachment
All parts (excluding the oxygen tubing & oxygen reservoir) are fully autoclavable up to 134°C
Transparency and resilience is better



LB3120
SEBS Resuscitator

Single patient use only.
Built-in 60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, PVC mask and oxygen tubing.
100% Medical level SEBS material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30mm exhale port for PEEP or filter attachment.



LB3130
PVC Resuscitator

Single patient use only.
Built-in 60cm H₂O pressure relief valve.
The set includes oxygen reservoir bag, PVC mask and oxygen tubing.
100% Medical level PVC material.
Latex-free components.
All-in-one intake valve connects to oxygen reservoir bag directly.
Additional accessories and customized label are available.
Built-in 30mm exhale port for PEEP or filter attachment.



LB3150
CPR Mask

Featurers:
Non-reusable, single use only.
Foldable or non-foldable mask desgins.
High transparency for better visualization.
100% medical level PVC material.
Latex-free components.

Specification:
Foldable cushion mask with oxygen port.
One way valve with high efficiency filter.
Elastic head trap.
Individual colorful hard case package.
With Individual PP box.
Box color: orange, blue, green, yellow, white.



LB3210
Disposable Silicone Laryngeal Mask

Medical grade silicone for super compatibility and patient comfort
are available with seven sizes for all patient use

Cat. NO.	Type	Applicable People
LB321010	1.0	< 5kg
LB321015	1.5	5-10kg
LB321020	2.0	10-20kg
LB321025	2.5	20-30kg
LB321030	3.0	30-50kg
LB321040	4.0	50-70kg
LB321050	5.0	> 70kg



LB3220
Reusable Silicone Laryngeal Mask

Reusable Silicone Laryngeal Mask are designed to be sterilized and used up to 40 times
Durable high temperature resistant material suitable for autoclave sterilization
Additional user card for record of device usage times and patient information is supplied in each pouch

Cat. NO.	Type	Applicable People
LB322010	1.0	< 5kg
LB322015	1.5	5-10kg
LB322020	2.0	10-20kg
LB322025	2.5	20-30kg
LB322030	3.0	30-50kg
LB322040	4.0	50-70kg
LB322050	5.0	> 70kg



LB3230
Disposable PVC Laryngeal Mask

Medical grade PVC with good compatibility
Available with seven sizes for all patient use

Cat. NO.	Type	Applicable People
LB323010	1.0	< 5kg
LB323015	1.5	5-10kg
LB323020	2.0	10-20kg
LB323025	2.5	20-30kg
LB323030	3.0	30-50kg
LB323040	4.0	50-70kg
LB323050	5.0	> 70kg

Oropharyngeal Airways

Comfortable and safe hard non-toxic polyethylene material with color-coded bite block for easy identification

LB3310 Guedel Airways

Cat. NO.	Size	Color Code
LB3310030	0000#, 30mm	Natural
LB3310040	000#, 40mm	Pink
LB3310050	00#, 50mm	Blue
LB3310060	0#, 60mm	Black
LB3310070	1#, 70mm	White
LB3310080	2#, 80mm	Green
LB3310090	3#, 90mm	Yellow
LB3310100	4#, 100mm	Red
LB3310110	5#, 110mm	Light Blue
LB3310120	6#, 120mm	Orange



LB3320 Berman Airways

Two side channel design provides passage of suction catheter to clear secretions and also provides space for oral care.

Packing: Packed in individual PE bag or blister packing or box.

Cat. NO.	Size
LB3330040	000#, 40mm
LB3330050	00#, 50 mm
LB3330060	0#, 60 mm
LB3330070	1#, 70 mm
LB3330080	2#, 80 mm
LB3330090	3#, 90 mm
LB3330100	4#, 100mm
LB3330110	5#, 110 mm



LB3410

Nasopharyngeal Airway

Made of non-toxic PVC, soft and comfortable
Only for nasal airway

Cat. NO.	Size
LB3410025	2.5 mm
LB3410030	3.0 mm
LB3410035	3.5 mm
LB3410040	4.0 mm
LB3410045	4.5 mm
LB3410050	5.0 mm
LB3410055	5.5 mm
LB3410060	6.0 mm
LB3410065	6.5 mm
LB3410070	7.0 mm
LB3410075	7.5 mm
LB3410080	8.0 mm
LB3410085	8.5 mm
LB3410090	9.0 mm
LB3410095	9.5 mm



Endotracheal Tube

X-ray opaque line throughout the length of the tube
Non-toxic PVC, transparent, soft and smooth

Standard Endotracheal Tube

LB5010 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5010020	2.0	LB5010070	7.0
LB5010025	2.5	LB5010075	7.5
LB5010030	3.0	LB5010080	8.0
LB5010035	3.5	LB5010085	8.5
LB5010040	4.0	LB5010090	9.0
LB5010045	4.5	LB5010095	9.5
LB5010050	5.0	LB5010100	10.0
LB5010055	5.5		
LB5010060	6.0		



LB5010C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5010025C	2.5	LB5010065C	6.5
LB5010030C	3.0	LB5010070C	7.0
LB5010035C	3.5	LB5010075C	7.5
LB5010040C	4.0	LB5010080C	8.0
LB5010045C	4.5	LB5010085C	8.5
LB5010050C	5.0	LB5010090C	9.0
LB5010055C	5.5	LB5010095C	9.5
LB5010060C	6.0	LB5010100C	10.0



Oral Preformed Endotracheal Tube

LB5020 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5020020	2.0	LB5020065	6.5
LB5020025	2.5	LB5020070	7.0
LB5020030	3.0	LB5020075	7.5
LB5020035	3.5	LB5020080	8.0
LB5020040	4.0	LB5020085	8.5
LB5020045	4.5	LB5020090	9.0
LB5020050	5.0	LB5020095	9.5
LB5020055	5.5	LB5020100	10.0
LB5020060	6.0		



LB5020C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5020030C	3.0	LB5020070C	7.0
LB5020035C	3.5	LB5020075C	7.5
LB5020040C	4.0	LB5020080C	8.0
LB5020045C	4.5	LB5020085C	8.5
LB5020050C	5.0	LB5020090C	9.0
LB5020055C	5.5	LB5020095C	9.5
LB5020060C	6.0	LB5020100C	10.0
LB5020065C	6.5		



Endotracheal Tube

X-ray opaque line throughout the length of the tube
Non-toxic PVC, transparent, soft and smooth

Nasal Preformed Endotracheal Tube

LB5030 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5030020	2.0	LB5030065	6.5
LB5030030	3.0	LB5030070	7.0
LB5030035	3.5	LB5030075	7.5
LB5030040	4.0	LB5030080	8.0
LB5030045	4.5	LB5030085	8.5
LB5030050	5.0	LB5030090	9.0
LB5030055	5.5	LB5030095	9.5
LB5030060	6.0	LB5030100	10.0

LB5030C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5030030C	3.0	LB5030065C	6.5
LB5030035C	3.5	LB5030070C	7.0
LB5030040C	4.0	LB5030075C	7.5
LB5030045C	4.5	LB5030080C	8.0
LB5030050C	5.0	LB5030085C	8.5
LB5030055C	5.5	LB5030090C	9.0
LB5030060C	6.0	LB5030095C	9.5
		LB5030100C	10.0

Reinforced Endotracheal Tube

LB5040 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5040030	3.0	LB5040065	6.5
LB5040035	3.5	LB5040070	7.0
LB5040040	4.0	LB5040075	7.5
LB5040045	4.5	LB5040080	8.0
LB5040050	5.0	LB5040085	8.5
LB5040055	5.5	LB5040090	9.0
LB5040060	6.0	LB5040095	9.5

LB5040C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5040030C	3.0	LB5040065C	6.5
LB5040035C	3.5	LB5040070C	7.0
LB5040040C	4.0	LB5040075C	7.5
LB5040045C	4.5	LB5040080C	8.0
LB5040050C	5.0	LB5040085C	8.5
LB5040055C	5.5	LB5040090C	9.0
LB5040060C	6.0	LB5040095C	9.5

Tracheostomy Tube

Made of non-toxic PVC
Used for tracheostomy

Tracheostomy Tube

LB5110 Uncuffed

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030	3.0	LB5110070	7.0
LB5110035	3.5	LB5110075	7.5
LB5110040	4.0	LB5110080	8.0
LB5110045	4.5	LB5110085	8.5
LB5110050	5.0	LB5110090	9.0
LB5110055	5.5	LB5110095	9.5
LB5110060	6.0	LB5110100	10.0
LB5110065	6.5		

LB5110C High Volume Low Pressure Cuff

Cat. NO.	Size ID (mm)	Cat. NO.	Size ID (mm)
LB5110030C	3.0	LB5110070C	7.0
LB5110035C	3.5	LB5110075C	7.5
LB5110040C	4.0	LB5110080C	8.0
LB5110045C	4.5	LB5110085C	8.5
LB5110050C	5.0	LB5110090C	9.0
LB5110055C	5.5	LB5110095C	9.5
LB5110060C	6.0	LB5110100C	10.0
LB5110065C	6.5		

SH0636 Endotracheal Intubation Set

Specification:

Endotracheal tube with stylet 1pc
Suction catheter FR.14 with finger controller tip 1pc
Suction catheter FR.16 with finger controller tip 1pc
Gloves 1pair
Guedel airway 1pc
Mouth pad 1pc
Bite block 1pc
Syringe 5cc
Gauze swab 2pcs
Suction connecting tube 1pc
Face mask 1pc
Cotton balls with silicone oil 2pcs

Endotracheal Tube Introducer (Bougie)

Smooth elbow design of tip can assist intubation
decreasing the damage of soft tissue
scale mark indicates distance to the tip
recommend for difficult intubation and replacing
an endotracheal tube
Available with solid type and hollow type.
Hollow type can facilitate delivery of oxygen in emergency


Cat. NO.	Size (Fr)	Type	Applicable ETT
LB521006	6	Solid	< 3.5mm
LB521010	10	Solid	4.0-5.0mm
LB521015	15	Solid	> 5.5mm
LB522010	10	Hollow	4.0-5.0mm
LB522015	15	Hollow	> 5.5mm



LB5310
Intubating Stylet

Made of malleable aluminum covered with plastic sheath
Soft rounded tip for minimizing patient trauma during intubation.
Flexible use for assistance of intubation.

Cat.No.	Size ID(Fr)	Appliable ETT	Length(mm)
LB531006	6	2.5-4.5	302
LB531010	10	4.0-6.0	382
LB531012	12	5.5-6.5	382
LB531014	14	6.0-10.0	382



LB6200
Suction Catheter

Made of medical grade PVC

Product Description :
Tube is frosted or clear
Open distal end with 2 side holes
With X-ray or without X-ray

Choice of connectors :
LB6210 Plain type
LB6220 Cap-cone
LB6230 Y-type
LB6240 T-type
LB6250 Y-Type
Specification: Transparent Fr5-Fr22
Standard Tube Length: 50cm
Different length is available



LB6310
Feeding Tube

Made of medical grade PVC

Product Description :
Tube is frosted or clear
Closed distal end with 2 side holes
With X-ray or without X-ray
Scale line

Specification: Fr4-Fr8
Standard Tube Length: 40/45cm
Specification: Fr10-Fr22
Standard Tube Length: 120cm



LB6410
Stomach Tube

Made of medical grade PVC or 100% silicone

Product Description :
Tube is frosted or clear
Open distal end with 4 side holes
With X-ray or without X-ray
Scale line

Specification: Fr5-Fr22
Standard Tube Length: 120cm or 125cm



LB6510
Rectal Tube

Made of medical grade PVC

Product Description :
Tube is frosted or clear
Closed distal end with 2 side holes
Specification: Fr18-Fr36
Standard Tube Length: 40cm



LB6610
Nelaton Catheter (PVC)

Made of medical grade PVC

Product Description :
Tube is Play frosted or clear
Closed distal end with 2 side holes

Specification: Fr6-Fr24
Standard Tube Length: **LB6610** Male 40cm
LB6620 Female 20cm



Available with different connectors

LB6122
Extractor Mucus

Volume: 25ml
With or without filter

Size ID(Fr)	Cat.No.:LB6122 (Plain Type)	Cat.No.:LB6112 (Cap-Cone)
Fr6	LB612206	LB611206
Fr8	LB612208	LB611208
Fr10	LB612210	LB611210
Fr12	LB612212	LB611212
Fr14	LB612214	LB611214



LB6810
Closed Suction Catheter

Sputum suction pipe sealing film produced by imported resin,
strong toughness, feel is exquisite and conducive to health care workers
Sputum suction tube body design scale mark line,
facilitate throughout the sputum suction effect
Closed suction phlegm tube adopts full closed type design, can
be used to sync with breathing machine, at the time of sputum
suction oxygen to sustainable, convenient and safe operation,
to avoid cross-infection
Specification : Fr6-Fr16
Type: 24 hours/72 hours



URINARY


Urethral Catheter

- Made of Latex, silicone coated
- Color-coded for size identification
- Sterile by EO
- For single use only
- Plastic & rubber valves are available



SH0101
2-Way Latex Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
SH010106	6	3	Light Red	27
SH010108	8	3-5	Black	27
SH010110	10	3-5	Grey	27
SH010112	12	5-10/30	White	40
SH010114	14	5-10/30	Green	40
SH010116	16	5-10/30	Orange	40
SH010118	18	5-10/30	Red	40
SH010120	20	5-10/30	Yellow	40
SH010122	22	5-10/30	Violet	40
SH010124	24	30	Blue	40
SH010126	26	30	Pink	40



SH0103
2-Way Female Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
SH010312	12	5-10/30	White	27
SH010314	14	5-10/30	Green	27
SH010316	16	5-10/30	Orange	27
SH010318	18	5-10/30	Red	27
SH010320	20	5-10/30	Yellow	27
SH010322	22	5-10/30	Violet	27



SH0104
2-Way Tiemann Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
SH010412	12	5-10/30	White	27
SH010414	14	5-10/30	Green	27
SH010416	16	5-10/30	Orange	27
SH010418	18	5-10/30	Red	27
SH010420	20	5-10/30	Yellow	27
SH010422	22	5-10/30	Violet	27



SH0102
3-Way Latex Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
SH010216	16	30	Orange	40
SH010218	18	30	Red	40
SH010220	20	30	Yellow	40
SH010222	22	30	Violet	40
SH010224	24	30	Blue	40
SH010226	26	30	Pink	40


always by your side

URINARY



SH0105
2-Way All Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
SH010506	6	3	Light Red	31
SH010508	8	3-5	Black	31
SH010510	10	3-5	Grey	31
SH010512	12	5-10/30	White	40
SH010514	14	5-10/30	Green	40
SH010516	16	5-10/30	Orange	40
SH010518	18	5-10/30	Red	40
SH010520	20	5-10/30	Yellow	40
SH010522	22	5-10/30	Violet	40
SH010524	24	30	Blue	40
SH010526	26	30	Pink	40




SH0106
3-Way Silicone Foley Catheter

Cat. NO.	Size (Fr/Ch)	Balloon Capacity (cc)	Color Code	Length(cm)
SH010616	16	30	Orange	40
SH010618	18	30	Red	40
SH010620	20	30	Yellow	40
SH010622	22	30	Violet	40
SH010624	24	30	Blue	40
SH010626	26	30	Pink	40



SH0107
Urethral Catheter (Red Latex)

Cat. NO.	Size (Fr/Ch)	Length(cm)
SH010706	6	27
SH010708	8	27
SH010710	10	27
SH010712	12	40
SH010714	14	40
SH010716	16	40
SH010718	18	40
SH010720	20	40
SH010722	22	40
SH010724	24	40



LB6200L
Latex Suction Tube

Cap-cone Suction Connector
Latex Tube
Standard with cap-cone connector
Available with different types of connectors
Total length: 40cm



SH0108
Malecot Catheter

Cat. NO.	Size (Fr/Ch)	Length(cm)
SH010812	12	40
SH010814	14	40
SH010816	16	40
SH010818	18	40
SH010820	20	40
SH010822	22	40
SH010824	24	40
SH010826	26	40
SH010828	28	40
SH010830	30	40



SH0110
Condom Catheter

Made of natural latex
By connecting pipe, fixture, membrane,
adhesive strap, penis
Size: 25mm, 30mm, 35mm, 40mm



SH0109
Latex T Tube Drainage

Made of natural latex, and silicone coated
Used for open gallbladder
The stem is often connected to a bile bag
Size: Fr10-Fr26



SH0110S
Silicone Condom Catheter

Made of 100% imported medical grade silicone
Pop On: suitable for old person or sports person.
Standard: suitable for regular length and adhesive for daily wear.
Wide Band: More adhesive and viscous force
than standard for added security.

Cat.No.	Type	Diameter				
		S	M	L	XL	XXL
SH0110SP	Pop On	24mm	28mm	31mm	35mm	40mm
SH0110SS	Standard	24mm	28mm	31mm	35mm	40mm
SH0110SW	Wide Band	24mm	28mm	31mm	35mm	40mm



SH0201
Economic Urine Bag

Capacity: 1000ml/1500ml/2000ml
without bottom outlet
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: individual PE bag, 10pcs/inner bag, 250pcs/ctn



SH0202
Economic Urine Bag

Capacity: 2000ml
Push Pull Valve
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: individual PE bag, 10pcs/inner bag, 250pcs/ctn



SH0205
Urinary Drainage Bag

Capacity: 2000ml
Twist Turn Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Sterilization: EO
Packing: individual PE bag, 100pcs/ctn




SH0206
Urinary Drainage Bag

Capacity: 2000ml
T- Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
Sterilization: EO
Packing: individual PE bag, 100pcs/ctn



SH0203
Economic Urine Bag

Capacity: 2000ml
Twist Turn Valve
With bottom outlet
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: individual PE bag, 10pcs/inner bag, 250pcs/ctn



SH0204
Economic Urine Bag

Capacity: 2000ml
T-Tap Valve
with bottom outlet
Adapter with cap
Anti-Reflux Flutter Valve
Inlet tube: O.D 6.5mm; 90cm Length
Sterilization: EO
Packing: individual PE bag, 10pcs/inner bag, 250pcs/ctn



SH0207
Urinary Drainage Bag

Capacity: 2000ml
T- Valve
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
With bed sheet clamp
With universal double-hook
Sterilization: EO
Packing: individual PE bag, 100pcs/ctn



SH0208
Urinary Drainage Bag

Capacity: 2000ml
With Single Hanger
Anti-Reflux Flutter Valve
Inlet tube: O.D 10mm; 100cm Length
Needleless Sampling Port with Cap
With Slide Clamp
Silicone Outlet Tube with Plastic Folding Clamp
Sterilization: EO
Packing: individual PE bag, 40pcs/ctn



SH0209
Urinary Drainage Bag

Capacity: 2000ml
Twist Turn Valve
With ribbon
Anti-Reflux Flutter Valve
Inlet tube: O.D 9mm; 100cm Length
Adapter with cap
Sterilization: EO
Packing: individual PE bag, 100pcs/ctn



SH0210
Urinary Drainage Bag

Capacity: 2000ml
with bed sheet clamp
Anti-Reflux Flutter Valve, Air-Vent on the Bag
Inlet tube: O.D 9mm; 100cm Length
Needleless Sampling Port with Cap
With Single Hanger and Rope Hanger
Silicone Outlet Tube with Metal Folding Clamp
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0213
Urinary Drainage Bag

Capacity: 2000ml
T-Tap Outlet Valve, Rigid Dog House
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Anti-Reflux Drip Chamber, Air-Vent on the Bag
Universal Double-Hook and Rope Hanger
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0214
Urinary Drainage Bag

Capacity: 2000ml/2500ml
T-Tap Outlet Valve, Rigid Dog House
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Three Parts Anti-Reflux Drip Chamber
air-Vent on the Bag
Universal Double-Hook
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0211
Urinary Drainage Bag

Capacity: 2000ml
T-Tap Outlet Valve, Rigid Dog House
With bed sheet clamp, with Slide Clamp
Inlet tube: O.D 10mm; 110cm Length
Needleless Sampling Port with Cap
Three Parts Anti-Reflux Drip Chamber
Air-Vent on the Bag
Universal Double-Hook and Rope Hanger
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0212
Urinary Drainage Bag

Capacity: 2000ml
T-Tap Outlet Valve
With plastic cover, with bed sheet clamp, with Slide Clamp
Anti-Reflux Drip Chamber, air-Vent on the Bag
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
Universal Double-Hook
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0215
Urinary Drainage Bag

Capacity: 4000ml/5000ml
T-Tap Outlet Valve, Rigid Dog House
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Two Parts Anti-Reflux Drip Chamber
air-Vent on the Bag
Universal Double-Hook and Rope Hanger
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0216
Urinary Drainage Bag

Capacity: 2000ml
T-Tap Outlet Valve
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Two Parts Anti-Reflux Drip Chamber
air-Vent on the Bag
Sterilization: EO
Packing: individual blister packing, 40pcs/ctn



SH0217
Urinary Drainage Urine Meter

400cc Two-Chamber Urine meter, 2000ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Universal Double-Hook and Rope Hanger
Silicone Outlet Tube with Metal Folding Clamp
and Splash Guard
Sterilization:EO
Packing: individual blister packing, 10pcs/ctn



SH0218
Urinary Drainage Urine Meter

400cc Two-Chamber Urine meter, 2500ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Universal Double-Hook
T-Tap Outlet Valve, with plastic cover
Sterilization: EO
Packing: individual blister packing, 10pcs/ctn



SH0219
Urinary Drainage Urine Meter

500cc Urinometer with Drip Chamber
2600ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Needleless Sampling Port with Cap
with bed sheet clamp, with Slide Clamp
Universal Double-Hook and Rope Hanger
T-Tap Outlet Valve
Sterilization: EO
Packing: Individual blister Packing, 10pcs/ctn



SH0220
Urinary Drainage Urine Meter

500cc Urine meter with 2000ml Urine Bag
Inlet tube: O.D 10mm; 120cm Length
Drip Chamber with Air-Vent and Anti-Reflux Valve
Needleless Sampling Port with Cap, Plastic Hanger
bed Sheet Clamp, with Slide Clamp
T-Tap Outlet Valve
Sterilization: EO
Packing: Individual blister Packing, 10pcs/ctn



SH0221
Urinary Drainage Leg Bag

500ml/750ml Leg Bag
Anti-Reflux Flutter Valve
Twist Valve
Two Pre-Assembled Straps
Sterilization: EO
Packing: individual PE Bag, 48pcs/ctn



SH0222
Urinary Drainage Leg Bag

500 ml, non-woven fabric on the back
Small and white T valve plug
Inlet tube: O.D 9mm; 30cm Length
Silicone tube 4.5 cm
Sterilization: EO
Packing: individual blister packing, 100 pcs/ctn




SH0223
Urinary Drainage Leg Bag

750 ml, non-woven fabric on the back, three cavity
Small and white T valve plug,
Inlet tube: O.D 9mm; 30cm Length
Needless sampling port with cap
Silicone tube 4.5cm
Two straps are available
Sterilization: EO
Packing: individual blister packing, 100 pcs/ctn



SH0224
Urinary Drainage Leg Bag

750 ml, non-woven fabric on the back
Small and white T valve plug,
Inlet tube: O.D 9mm; 30cm Length
Needless sampling port with cap
Two Pre-assembled Straps
Sterilization: EO
Packing: individual blister packing, 100 pcs/ctn



SH0229
Pediatric Urinary Collection Bag

Safety Adhesive, PVC Bag, Sponge
Capacity: 100ml, 200ml
Sterilization: EO
Packing: individual PE bag, 100 pcs/inner box, 5000pcs/ctn



SH0244
Closed Wound Drainage System Kit

Made of imported medical grade silicone
This product includes 100cc/150cc/200cc/400cc Silicone Reservoir 7mm/10mm Flat drains or other round drains and stainless steel trocar
Can use separately or using Trocar for insertion and combining with silicone Reservoir for fluids suction and collection
Used for various surgeries
Latex-free
Sterilization: EO
Packing: individual blister packing, 20pcs/ctn




SH0242 Closed Wound Drainage System (spring)
SH0243 Closed Wound Drainage System (hollow)

Our closed wound drainage system, including 400cc, 3-spring evacuator, PVC tubing and "Y"connector, and drain with trocar.The 3-spring evacuator design is minimize collapse-resistance, packed in double liner Tyvek, to make sure medical cleanness on O.R.operation.
Capacity: SH0242 (200ml, 400ml, 800ml)
SH0243 (400ml, 500ml, 800ml)
Silicone or PVC tube, 120cm Length
Sterilization: EO Packing: Individual blister packing, 20pcs/ctn





SH0234
Enteral Feeding Bag Set

Volume: 500ml, 800ml,1000ml, 1200ml, 1500ml
Gravity Type/Pump Type
Enteral Feeding Bag has rigid neck for easy filling and handling.
Strong, dependable hanging ring.
Easy-to-read graduation and easy-view translucent bag.
The exit port in the bottom allows completely drainage.
Sterilization: EO
Packing: individual PE bag, 24pcs/ctn



SH0236
EVA TPN (total parenteral nutrition) Bag

Features:
- Non-PVC material, DEHP free, can be used to fat emulsion and other fat soluble liquid and medicine.
- No adverse effects on infants and young children and pregnant women.
- Flexible tube and bag to avoid kinking and breaking.

Item	Volume	Details
	500ml	Eva bag body Eva inlet tube
	1000ml	
	2000ml	
	3000ml	
	500ml	Eva bag body
	1000ml	
	2000ml	
	3000ml	



SH0245
Disposable Drainage Bottle

Single bottle, 1600ml
Double bottle, 2000ml
Treble bottles, 2500ml

Constitution:

By attachment, drive pipe, bottle cap, bottle body, dip pipe, support and suspension hook constitution.

SH0241
Suction Liner



Product Characteristic:
Simplicity
Safety
Economical

Indication for Use:
This product is used for collecting the biohazardous fluids which produced in ICU and operation room.

Product Features:
1.Double safety design to ensure that no liquid is suctioned back to the tubing and suction unit, effectively prevent cross infection and protect the medical staff safe.
2.It's made of environment-friendly polyethylene. It's safe and non-toxic.
3.Easy to use, can be a single use, also can use with the stander.
4.Available in some specifications: with or without shutoff filter, with or without solidifier, 1000ml, 1500ml, 1800, 2000ml, 3000ml, 3500ml.

Product Advantages:
1.Soft and transparent(Convenient to observe liquid waste)
2.Double safety design(To prevent fluid into the suction systems by mistake)
3.It can be connected with relevant suction device or vacuum device, soft bag is seal and reliable inside and can be withstand the weight of 60kg when fully filled with liquid(Prevent infections liquid overflow caused by medical personnel)



Specification

Item	Description
Stand	Four tandem stand 4 canisters (a set)
	Two tandem stand 2 canisters (a set)
Hanger	Two tandem hanger stand
	Single hanger stand
	Single canister

SH030
New Style Ostomy Bag

Our ostomy bags are of excellent quality and advanced usage,making the wound healed in the good condition
Ostomy bag applications:collect excrement and urine after colostomy, ileostomy and urostomy.

Features:
1,The hydrocolloid skin barrier has the performance of anti-maceration, low allergy,and warping and leakage.
It has excellent exduation-absorbent and long sticking ability(up to 7 days).And it is soft and comfortable.
2,Hydrocolloid skin barrier combines with the highly permeable PU film makes it stick for a long time, no leakage of exudates,convenient to use.
3,The material of ostomy bag possesses high obstructing ability,it is comfortable and covert.
4,The one-piece system design discharge interstinal gas easily and filter effectively it excellent.



SH0301 Simple Type
SH0302 Foam Type
SH0303 Glue Pan Type

SH0304 Two-piece Type
SH0305 Foam Type
SH0306 Two-piece Type

HYPODERMIC



SH0601 Disposable Hypodermic Syringe

1ml, 2/3ml, 5ml, 10ml, 20ml, 30ml, 50/60ml
Three-part syringes/Two-part syringes (luer slip,
special design to prevent plunger slipping out)
Luer slip/luer lock
Latex/latex free piston
Sterilization: EO
PE/Blister individual package, 100pcs/Inner Box
W/needle, W/O needle



SH0608 Disposable Hypodermic Needle

Size: 18-30G
Length of needle: 1/2", 5/8", 3/4", 1" , etc.
Blister individual package
Sterilization: EO



SH0602 Auto Self-destruction Syringe

0.5mm, 0.1ml, 0.3ml, 0.5ml, 1ml
Fixed needle 27G-30G
Sterilization: EO
PE/Blister individual package, 200pcs/Inner Box



SH0604 Safety Clip Syringe

1ml, 2/3ml, 5ml, 10ml, 20ml
Safety syringe with protective cap
Luer lock
After injection, push the clip over the needle and lock it
onto the clip base
Sterilization: EO



SH0605
Tubercle Bacillus Syringe

1ml
luer slip
Latex/latex free piston
PE/Blister individual package
Sterilization: EO




SH0606
Insulin Syringe

0.3ml, 0.5ml, 1ml
Fixed needle 27G-31G
PE/Blister individual package
Sterilization: EO



SH0501
Disposable Infusion Set

For infusion bags and bottles
15/20/60 drops = 1 ± 0.1 ml
Small Chamber/Big Chamber
Latex tube
With Filter
Different types of components and parts can be changed as per customer's requirements
For gravity use
EO sterile
PE/Blister individual package



SH0502
Disposable Infusion Set

For infusion bags and bottles
15/20/60 drops = 1 ± 0.1 ml
Big Chamber
Latex free Y Site
With Filter
With micro-flow regulator
Different types of components and parts can be changed as per customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package



SH0401
Irrigation Syringe

50ml, 60ml, 80ml
Flat/Bulb/Ring
PE/Blister individual package
Sterilization: EO



SH0641
Enteral Syringe

The Oral/Enteral Dispenser are assembled by barrel, plunger, piston.
All parts and material for this product meet medical requirements, after sterilized by ETO.
The Oral/Enteral dispenser is used for delivering drug or food to oral or enteral.
Enteral Syringes are available in 8 sizes - 1, 3, 5, 6, 10, 20, 50 and 60ml.



SH0505
Disposable Blood Transfusion Set

For blood bags
Different types of components and parts can be changed as per customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package



SH0504
Burette Infusion Set

Burette Capacity: 100ml, 150ml
For infusion bags and bottles
60 drops = 1 ± 0.1 ml
with Fitter
with Y Site
Different types of components and parts can be changed as per customer's requirements
For gravity use
EO gas sterile
PE/Blister individual package



SH0607
Scalp Vein Set

Size: 18-30G
Length of needle: 1/2", 5/8", 3/4", 1", etc.
Tube length: 30cm
Luer slip/luer lock
PE/Blister individual package
Sterilization: EO



SH0610
Insulin Pen Needle

Ultra fine cannula for maximum patient comfort, patient point with triple beveling for easy and pain decreases
The needle is fit for Humapen (EliLily, USA)
NovopenIII (Novo Nordisk, Denmark)
Sterile, non-pyrogenic
Sterilization: EO



SH050901
Needle Free Connector

Features:
1.Transparent housing material: Polycarbonate or Copolyester.
2. Metal-free & MRI-compatible.
3. Latex-free.
4. Multiple insertion at least 100 times/day.
5. Priming Volume: 0.09ml.
6. Perfect flowing rate: 350 ml/min under one meter water pressure, tested by Baihe Technical Team.



SH050902
Y-port Valves

Material: Polycarbonate or Copolyester
Size: ID 4.01mm, 3.6mm



SH0701
Disposable Blood Bag
For the collection of a blood volume: 250, 350, 450, 500ml.
Is filld with 63ml of CPDA-1 ANTICOAGULANT SOLUTION
Each 100ml CPDA-1(U.S.P.) contains:
Citric Acid (anhydrous)-----0.299g
Sodium Citrate (dihydrate)-----2.63g
Monobasic Sodium Phosphate (monohydrate)-----0.222g
Dextrose (monohydrate)-----3.19g
Adenine-----0.0275g
Water for injection, a sufficient quantity to make-----100ml
Needle: 16 gauge
Steam sterilization
Sterile, non-toxic, non-pyrogenic fluid path



SH0507



SH0508

SH0507
Heparin Cap
Material: PC
Sterilization: EO
Individal Hard Blister Packing, 100pcs/box, 5000pcs/ctn

SH0508
Medical Disposable Consumable Spiral Stoppers
Material: ABS Sterilization: EO
Color: Red, White, Blue, etc.
Individual Hard Blister Packing, 100pcs/box



SH050903
Tube End Valves
Material: Polycarbonate or Copolyester
Size: ID 2.0mm, 2.5mm, 2.8mm, 3.7mm, 4.1mm



SH050904
Positive Pressure Needle Free Connector
Features:
1. Patent design of positive pressure avoids back flow of blood during removal of syringe. This will help to prevent blood clot at the catheter tip inside the blood vessel.
2. The housing is impregnated with PC Ag+ which helps to lower infection.
3. Protruding design of the stalk at the injection port follows infection control protocol.
4. High quality spring ensures multiple insertion at the injection site without leaking.
5. Two sealing rings at the upper and lower end of the stalk keeps the connector out air, fluid and outside substance.
6. The straight flow of fluid path creates less turbulence and complies with proper transfusion protocol.
7. The holding capacity of remaining fluid is less than 0.06ml.



SH0506
Three Way Stopcock

Parts: Body(PC), Plas(PE), Screw Cap(ABS)
Individual Hard Blister Packing, 50pcs/box, 1000pcs/ctn
Sterilization: EO



SH051001
Three way stopcock with Extension Tube

The length of extension tube can be changed as customer's requirements
Different types of components and parts can be changed as customer's requirements
Sterilization: EO
Packing: Individual blister Packing
Tube Length:
10-25cm, 50pcs/ctn
30-50cm, 50pcs/ctn
80cm, 50pcs/ctn
100-120cm, 40pcs/ctn



SH0511
Needle Free Adapters For Infusion Set

- 1.This device is designed to meet the needs of general I.V. Therapy, Anaesthesia Cardiovascular,
- 2.ICU&CCU, Recovery&Oncology.
- 3.Safe for infusion IV medication including most high risk anti-neoplastics, or aspirating blood without risk of hemolysis.

- Catheter Diameter:**
- 1.One-Way Extension Tube with One Normal Needle Free Connector
 - 2.One-Way Extension Tube with One Positive Pressure Needle Free Connector
 - 3.Two-Way Extension Tube with Two Normal Needle Free Connectors
 - 4.Two-Way Extension Tube with Two Positive Pressure Needle Free Connectors
 - 5.Three-Way Extension Tube with Three Normal Needle Free Connectors
 - 6.Three-Way Extension Tube with Three Positive Pressure Needle Free Connectors
 - 7.Four-Way Extension Tube with Four Normal Needle Free Connectors
 - 8.Four-Way Extension Tube with Four Positive Pressure Needle Free Connectors



SH0511002

Screw Extension Tube



SH0511003

Mini-Quantity Infusion
Extension Tube (OD:2.0mm)



SH0511004

T-Extension Tube



SH0513
Micro-flow Regulator

Flow range: 5ml/h-250ml/h
Materials:
Regulator: ABS
Tube: PVC
Other accessories: PC & ABS
Extension Tube Length: 45cm is common, others are also available.
Sterilization: EO
Packing Details: Poly-bag per pcs, 200pcs/inner box, 20boxes/carton



SHA01
Central Venous Catheter

Advanced automatic drilling technology and equipment, a more accurate location and the pore diameter of the pipe, surface smooth and without cutting residues and injury.

Special custom polyurethane fine materials, delicate and smooth tube, welding head end, ultra smooth surface treatment and X line to show the position, make a catheter inserted into the blood vessels in more easy and safe.



SHA07
Hemodialysis Blood Tubing Set

Hemodialysis Blood Tubing Set Consists of arterial and venous lines, used during dialysis and are attached with A. V. fistula needle and dialyzer.

Features:

- 1. Medical grade materials displaying good biocompatibility.
- 2. Fluid pathway is smooth inside, no clotting or air bubbles.
- 3. Tubes are flexible, kink-resistant.
- 4. Color Coding.
- 5. Pump segment has excellent anti-fatigue performance.
- 6. Luer-lock connectors according to international standard.
- 7. Various series, good applicability.
- 8. Individual and customized services provided.



SHA03
PICC Kit

- 1. The connector is round and smooth with PASV valve.
A) Valve opens inward for sampling under negative pressure.
B) Valve opens outward for infusion under positive pressure.
C) Valve remains closed to prevent air embolism, blood reflux and clot under normal pressure
- 2. The catheter is made from premium silicone, which has excellent biocompatibility.
- 3. Adjustable length of catheter.
- 4. Detectable under X ray, the end tip enhances radiopacity.
- 5. The style of catheter ensures placement more convenient and smooth.
- 6. Variety of kit configuration, initial guide wire and positive, pressure needless connector is available for choice.



SHA05
Hemodialysis Catheter Kit

- Standard Kits include:**
- 1. Hemodialysis Catheter
 - 2. Vessel Dilator
 - 3. Introducer Needle
 - 4. Syringe
 - 5. Guide-wire
 - 6. Adhesive Wound Dressings
 - 7. Herparin Caps
 - 8. Scalpel
 - 9. Needle with Suture
- Optional Compound Kits include:**
- 1. Hemodialysis Catheter
 - 2. 5ml Syringe
 - 3. Surgical Gloves
 - 4. Surgical Pledget
 - 5. Surgery Sheet
 - 6. Surgery Towel
 - 7. Sterile Brush
 - 8. Gauze Pad



SHA06
Long Term Hemodialysis Catheter Kit

- Characteristics:**
- 1. The Catheter is made from premium silicone, which has excellent biocompatibility.
 - 2. The dacron cuff can protect immigration of catheter and reduce the chance of CRBI (Catheter Related Bloodstream Infections).
 - 3. The silicone material makes catheter more floppy, hence it can minimize the trauma to the vein and reduce the incident of clot.
- Standard Kits include:**
- 1. Silicone Catheter
 - 2. Vessel Dilator
 - 3. Valved Peelable Introducer
 - 4. Introducer Needle
 - 5. Guide-wire
 - 6. Tunneler (Tocar)
 - 7. Adhesive Wound Dessings
 - 8. Herparin Caps
 - 9. Gauze Sponges
 - 10. Scalpel
 - 11. Syringe
 - 12. Needle with Suture



SHA02
Central Venous Catheter Kit

- Standard Kits include:**
- 1. Central Venous Catheter
 - 2. Guide-wire
 - 3. Vessel Dilator
 - 4. Clamp
 - 5. Fastener: Catheter Clamp
 - 6. Introducer Needle
 - 7. Introducer Syringe
 - 8. Injection Needle
 - 9. Injection Cap
- Optional Compound Kits include:**
- 1. Central Venous Catheter
 - 2. Standard Kit Accessories
 - 3. Surgical Gloves
 - 4. Surgical Pledget
 - 5. Surgery Sheet
 - 6. Surgery Towel
 - 7. Sterile Brush
 - 8. Gauze Pad
 - 9. Suture of Needle
 - 10. Wound Dressing
 - 11. Scalpel



SHA04
Drainage Catheter Kit

- Standard Kits include:**
- 1. Drainage Catheter
 - 2. Guide-wire
 - 3. Vessel Dilator
 - 4. Introducer Needle
 - 5. Introducer Syringe
 - 6. Injection Needle
 - 7. Injection Cap
 - 8. Connecting tube
- Optional Compound Kits include:**
- 1. Drainage Catheter
 - 2. Standard Kit Accessories
 - 3. Surgical Gloves
 - 4. Surgical Pledget
 - 5. Surgery Sheet
 - 6. Surgery Towel
 - 7. Sterile Brush
 - 8. Gauze Pad
 - 9. Suture of Needle
 - 10. Wound Dressing
 - 11. Scalpel



SH0512
Disposable Pressure Transducer

Disposable Pressure Transducer is made of medical PC and PVC material.

Operating Parameters:

- Operating Pressure Range: -30mmHg~+300mmHg
- Operating Temperature Range: 15℃~40℃
- Storage Range: -25℃~+70℃
- Atmosphere Pressure Range: 56.7kPa~113kPa
- Zero Output: ±25mmHg
- Sensitivity: 5 μV/V/mmHg±2%
- Input Impedance: 1200 Ω~3200 Ω
- Output Impedance: 285 Ω~315 Ω
- Zero Thermal Drift: ±0.3mmHg/℃
- Zero Drift: ≤2mmHg
- Natural Frequency: 40Hz nominal for a standard kit >200Hz for transducer alone
- Sensitivity Thermal Drift: ±0.1%/℃
- Leak Current: <2 μA at 230VRMS, 50Hz

Note: rated voltage 6.00VDC; ambient temperature 25 ℃.



- Specification:**
- 130cm pressure tubing with M/F connector 1piece
 - 30cm pressure tubing with M/F connector 1piece
 - Disopsable pressure transducer 1piece
 - Flush device 3 cc/h (30cc/h) 1piece
 - Male non-vented Cap 2piecd
 - 150cm single IV set 1piece
 - 3-way stopcock 2piece



SH060901
Dental Needle

Size: 27g X 38mm, 30g X 38mm , etc.
Sterilization: EO
Packing: 100pcs/Box, 5000pcs/ctn



SH060902
Dental Needle

Size: 18G-25G
Packing: 100pcs/Bag,
600bags/ctn



SH0620 Spinal Needle

Size: 18g~28g x 90mm
Sterilization: EO
Packing: 50pcs/Box, 1000pcs/ctn




SH0617 A.V. Needle

Size: 15Ga, 16Ga, 17Ga.
Sterilization: EO
Packing: 1pc/bag, 10pcs/Box, 500pcs/ctn



SH0622 Epidural Set

Features:
Epidural needle: 16G/18G x 80/90mm
Epidural catheter: 1.0/0.8mm,clear marked, length longer than 850mm with 3 holes or more holes, closed tip
Epidural filter: 0.22um, Luer lock
LOR syringe: 7ml, 10ml available
Catheter adaptor: TBA, secure the catheter completely
Single use to avoid the cross infection
Use in the clinic of epidural nerve or subarachnoid anesthesia



SH0616
Safety Disposable I. V. Catheter

Cat.NO.	Long	Specification	Color
SH0614	45mm	18G	Green
	32mm	20G	Pink
	25mm	22G	Blue
	20mm	24G	Yellow
	15mm	26G	Purple



SH0615



SH0613



SH0614

I. V. Cannula
Normal/Safety
Pen type/With wings/With injection port
Size: 14G-26G
PTFE/FEP/PU flex catheter
Sterilization: EO
Packing: Hard Blister Packing

All the popular sizes are available in both a carbon or stainless steel sterile
Blade shapes #9 to 17 are used in histology units and are all compatible with the associated
#3handles (Blade shape 11 can be also used for gardening with a plastic or metal handle).
Blade shapes #18 to 36, #60 and #70 are used for major invasive procedures and are all
compatible with the associated #4 handles
Sterilization: Gamma
Packing: 1pcs/Sealed Aluminum Foil Bag,
100pcs/box, 5000pcs/ctn



SH0901
Disposable Surgical Blade

Disposable Electrosurgical Pencil (Tips)



SH1101 Disposable Electrosurgical Pencil



SH110101
Electrode Shape: Blade
Handle Length: 158mm
Hand control
Cable Length: 3m
Exposed electrode length: 40mm



SH110102
Electrode shape: Loop
Handle Length: 158mm
Cable Length: 3m
Hand control
suitable for gynecological operation
Electrode Shape: Loop
Electrode Exposed Length: 40mm



SH110103
Electrode Shape: Blade
Handle Length: 158mm
Foot control
Cable Length: 3m
Exposed electrode length: 40mm



SH0902 Stitch Cutters

All sizes are available in both a carbon or stainless steel sterile
Designed with a rounded bull nose tip, the Index family of Stitch
Cutter Blades are ideal for
cutting the wide range of modern sutures available today whilst
minimising patient discomfort
The range consists of a short standard stitch cutter with a larger
grip. The long handle enables
a more aseptic technique allowing good visibility whilst working
further from the wound site
Size: 65mm, 110mm
Sterilization: Gamma
Packing: 1pc/Sealed Aluminum Foil Bag,
100pcs/box, 5000pcs/ctn



SH0903 Surgical Handle

Index stainless steel metal handles are supplied
non-sterile and individually wrapped
Sterilization: Gamma
Packing: 1pc/Sealed Aluminum Foil Bag,
10pcs/box, 500pcs/ctn



SH0904 Surgical Scalpels

Standard sheath guarded and blade protected sterile
disposable scalpels offer convenience in a number of
user environments: on the ward, in emergency or trauma
units or even out in the field with paramedics
Material: Carbon Steel or Stainless Steel
Sterilization: Gamma
Packing: 1pc/sealed Aluminum Foil Bag,
10pcs/Box, 500pcs/ctn



SH0905 Steel Lancet

Low cost: simple, flat, steel lancet
Sample volume: excellent blood flow
High quality: hardened surgical steel
Minimum discomfort: modern asymmetric design
Sterile: sterilized by gamma radiation
Size: 41x5.5mm
Packing: Individual pack, 5pcs of packaged,
200pcs/Box, 20000pcs/ctn



SH0906 Twist Lancets

Conventional twist-off cap design
Sterilized by gamma-radiation
Smooth tri-bevel needle tip for comfortable sampling experience
Versatile with most lancing devices
Size: 21G, 23G, 26G, 28G, 30G
Packing: 100pcs/box, 20000pcs/ctn



SH0907 Twist Lancets (Exclusive to Roche Diagnostics)

5 twist lancets (28G) are used with Roche Accu-Chek Softclix Lancing Device and the needle tip is siliconized for easy and smooth motion during lancing
Size: 28G, 30G
Packing: 100pcs/box, 20000pcs/ctn



SH0912 Button Activated Safety Lancets

Needle tip is sealed with cap through injection moulding
Pre-loaded, button operated design
Single-use only, sterilised by gamma irradiation, cannot be reused
High quality, ultra sharp tri-bevel needle, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G



SH0908 Lancing Device (Releasable) SH0909 Lancing Device

The user can install the lancet with the device easily. And the device will work with the lancet smoothly without any blockage. The lancet won't fall on its own when to fire the device. Retracting automatically after firing, the launcher is safe and reliable.



SH0910 Safety Lancet

Ergonomic design offers a more secure hold when to activate the device
Tri-bevel needle tip ensures the sharpness and accuracy
The spring buffer the puncture making the procedure a more comfortable experience
The complete range of needle sizes from 18G to 32G ensure the accurate
Blood samples for all popular meters
The needle retracts immediately after use
Sterilization : Gamma
Size: 21G, 23G, 26G



SH0911 Pressure Activated Safety Lancets

Needle is fully shielded before and after use
Pre-loaded, pressure activated design
Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization : Gamma
Size: 21G, 23G, 26G



SH2501 Thermal Accident Blanket

Material: PET Film
Waterproof, use to offset hypothermic reactions, reflective surface provides high visibility, reflects sun and heat for cooling in hot climates, space-age design retains up to 90% of body heat, Assists in prevention of post-accident shock and first aid application
Size: 132cm×210cm, 140cm×220cm, 160cm×210cm or other size per customer demand
Color: Gold, Silver, Gold+Silver



SH2502 Splints

Aluminum plate Macromolecule foam material
Nontoxic and nonbenzenoid
Free of effect of temperature and climate
Type available: Fold, Unfold and Roll
Normal size: 5cm×12cm, 5cm×15cm, 7cm×11cm, 11cm×44cm, 11cm×91.5cm, 11cm×100cm



SH0921 / SH0923 / SH0924 Medical Razor

- Stainless steel blades • Single or double blade • All razors come with the highest quality Comfort-coated for smooth shaves



SH1001

Surgical Suture

Cat. NO.	Type
SH1001	Chromic Catgut
SH1002	Plain Catgut
SH1003	Silk Braided
SH1004	Nylon Monofilament
SH1005	Polyester Braided
SH1006	Polypropylene
SH1007	Polyglycolic Acid
SH1008	Polydioxanone (PDS/PDO)

Thread Length: 45cm, 60cm, 75cm, 90cm, 150cm etc.
 Thread Diameter: 8/0, 7/0, 6/0, 5/0, 4/0, 3/0, 2/0, 1/0, 1, 2, 3 etc.
 Needle Length: 6mm, 8mm, 12mm, 18mm, 35mm, 40mm, 50mm etc.
 Needle Curvature: Straight, 1/2circle, 1/2circle(double), 1/4circle, 1/4circle(double), 3/8circle, 3/8circle(double), 5/8circle, loop round
 Cross Section: round bodied, round bodied(heavy), curved cutting, curved cutting(heavy), reverse cutting, reverse cutting(heavy) etc.

SeaSight®

EXAMINATION



SH0802 Surgical Gloves

Size: 6.0#, 6.5#, 7.0#, 7.5#, 8.0#, 8.5#, 9.0#, 9.5#
Powdered/Powder free
1pair/pouch, gamma/EO sterile
AQL: 1.5
Packing: 50pairs/box, 500 pairs/ctn



SH0803 Latex Examination Gloves

Size: S,M,L
Powdered/Powder free
100pcs/box, non-sterile
AQL: 2.5
Packing: 100pcs/box, 1000 pcs/ctn



SH0804 Nitrile Examination Gloves

Size: XS, S, M, L
Powder free
100pcs/box, non-sterile
AQL: 1.5/2.5/4.0
Packing: 100pcs/box, 1000 pcs/ctn



SH0805 Vinyl Gloves

Size: S (3.0g,3.5g,4.0g) ; M (3.5g, 4.0g, 4.5g) ;
L (4.0g,4.5g,5.0g)
Powdered/Powder free
100pcs/box, non-sterile
AQL: 4.0
Packing: 100pcs/box, 1000 pcs/ctn
Corlor: Blue, White



SH0806 PE Gloves

Size: S,M,L
Thickness: 0.018mm, 0.02mm, 0.022mm etc.
Type: HDPE/LDPE
Packing: 100pcs/box, 10000 pcs/ctn



SH0807 Two Fingers Cots

Size: S,M,L
Packing: 100pcs/inner bag, 5000 pcs/ctn

Vaginal Speculum

always by your side



SH0811 Side Screw Type

Size: S,M,L



SH0812 Fastener Type

Size: S,M,L



SH0813 Spanish Type

Size: S,M,L



SH0814 Middle Screw Type

Size: S,M,L



SH0815 Push Pull Type

Size: S,M,L



SH0816 American Type

Size: S,M,L



SH0817 Disposable Anoscopes

Disposable anoscope with integrated LED light source
Single use, cable-less and self-illuminating proctoscope
designed to optimize mobility during diagnosis and treatment.
Button Battery SR-LR41, 1,5V
Battery life: 3 years
Light work time: 48 hours
Packing: individually sterilized packing



SH082101 Cervical Brush




SH082102 Cervical Brush



SH082103 Cervical Brush



SH082104 Cervical Brush



SH082105 Cervical Brush



SH082106 Cervical Scraper



SH082107 Cervical Smear



SH082108 Cervical Smear



SH0830
Wooden Tongue Depressor

Size:Adult (150mm × 18mm × 1.6mm)
Child (140mm × 14mm × 1.6mm)
EO Sterile / Non-Sterile
Packing: Single medical paper wrapped/in Bulk
100pcs/Box, 5000pcs/Ctn



SH0832
Cervical Scraper (Wooden)

EO Sterile / Non-Sterile
Packing : Single medical paper wrapped/in Bulk
100pcs/Box, 5000pcs/Ctn



SH0822
Gynecological Kit

Configuration :
1pc speculum
1pc Gynaecology pad(40*60cm)
1pair latex glove
1pc cervical smear



SH1201 Condom
Specification:
Surface: smooth, dotted, ribbed lubricated
Size: L, M, S (55mm±2、52mm±2、49mm±2).
Be available with spermacide nonoxynol 9
Be available with various fruit smell, such as apple/ banana/ strawberry/ cherry etc.
Be available with square rectangle circular, or butter packing.
If quantity is large enough, the condom is available with various Colors: such as natural/ red/pink/green etc.
Specification is available with En600,ISO4074, ASTM, WHO, SABS, etc.



SH1301
HCG Test

SH1301



Size: 2.5mm, 3.0mm, 3.5mm, 4.0mm, 4.5mm, 5.0mm, etc
Response time: 5 minutes
Sensitivity: 25mIU/ml
Specificity: 99%
Packing: 1pc/aluminum foil bag



SH1302 Urine Reagent Strips

The Urine Reagent Strips (URS) for Urinalysis are firm plastic strips to which several different reagent areas are affixed. Depending on the product being used, Urine Reagent Strips provide tests in urine for up to 14 different items. Items: Leukocytes;Nitrite;Urobilinogen;Protein;pH; Blood;Specific Gravity;Ketone;Bilirubin;Glucose; Ascorbate;Micro Albumin;Creatinine;Calcium.



SH1303 Dengue

The Dengue Rapid Test detection of IgG, IgM antibodies and/or NS1 antigen to dengue virus in human whole blood, serum or plasma as an aid in the diagnosis of dengue infections. Materials Provided: Dengue Rapid Test Cassettes Droppers Buffer Package insert



SH1304 Malaria

The Malaria Rapid Test use a rapid lateral flow chromatographic immunoassay in human blood specimen as an aid in the diagnosis of Malaria infection. Materials Provided: Malaria Rapid Test Cassettes Droppers Buffer Package insert



SH1305 Saliva Drug Test

The Oral Fluid Drug Test is a lateral flow chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in oral fluid. Materials Provided: One Step Multi-Drug Test Cubes/Cassettes Package insert Saliva collectors Security seal labels (for cubes use only) Color chart card for alcohol (when applicable)



SH1306 Urine Multi-drug Test

Specimen validity combinations available for dip card g cups: Creatinine (CRE), Glutaraldehyde (GLUT), Nitrite (NIT), Oxidants (OX), pH & Specific Gravity (SG). Materials Provided: One Step Multi-Drug Test Dip Cards/Cassettes/Cups Package insert Droppers (for cassettes use only) Disposable gloves (for cups use only) Procedure Card (for cups use only) Color chart cards for adulterant interpretation (when applicable) Color chart card for alcohol (when applicable)

SeaSight®

WOUND DRESSING



SH2001 Gauze Swabs

The pads are folded all by machine, pure 100% cotton yarn ensure the product soft and adherent. Superior absorbency makes the pads perfect for absorbing blood and exudates. In accordance with the customers' requirements we can produce different kinds of pads, such as folded, and folded, with x-ray and non x-ray. The adherent pads are perfect for operation. General specification: 2"×2", 3"×3", 4"×4" PLY: 8ply, 12ply, 16ply Mesh: 19×10, 19×15, 20×12, 26×18 EO Sterile / Non-Sterile



SH2002 Non-woven Swabs

Non-woven swabs from 70% polyethylene and 30% terylene. The pad is designed to wick away fluids and disperse them evenly. Product have been cut like "O" and "Y", so it is easy to use. It is mainly used to absorbing blood and exudates during operation and cleaning the wounds. General specification: 30g/m², 40g/m²; 2"×2", 3"×3", 4"×4" PLY: 4ply, 8ply EO Sterile / Non-Sterile



SH2005 Gauze Bandage

Material : 100% cotton
Cotton Yarn: 40's×40's
Weaving according to customer requirements, high pressure, high temperature degreasing, bleaching, refined but become
Packing : 1 roll/white paper, 1 doz blue kraft paper packing
General specification:
2" , 3" , 4" , 6" ×4.5m/5m/10m;
Mesh: 19×15, 20×12, 26×18
EO Sterile / Non-Sterile



SH2007 Cotton Roll

Absorbent cotton is pure 100% cotton fabric, degreased and bleached by advanced way. The cotton is perfect for beauty; light drainage wounds, the alcohol sterile, cleaning the medical instruments. Packing become rolls by white and blue paper or PE bag. General specification: 25g, 50g, 100g, 200g, 250g, 400g, 454g, 500g, 1000g EO Sterile / Non-Sterile



SH2003 Prewashed Laparotomy Sponge

Laparotomy Sponge for absorbing blood and exudates or for cleaning wound. Cotton Yarn: 40's x 40's It can prevent the wound from being infected. Our laparotomy sponges are sewn and cross-stitched. General specification: 45×45cm/30×30cm, PLY: 8ply, 16ply Mesh: 19×15/26×18 EO Sterile / Non-Sterile Prewashed/ Non-Prewashed



SH2004 Gauze Roll

Material : 100% cotton
Cotton Yarn: 40's×40's
Weaving according to customer requirements, high pressure, high temperature degreasing, bleaching, refined but become, with x-ray or without x-ray
Products without fluorescence, soft, has the very strong absorption ability, is the first choice for medical units use products, quality in line with the British pharmacopoeia standards, are widely used in other markets
General specification: 36"×100yard-4ply, 36"×100m-4ply, 36"×1000m-4ply; Mesh: 19×11/19×15



SH2006 Gauze Ball

100% pure cotton yarn of our product ensures a superior absorbency. There are two types of this product: with x-ray and non x-ray ball. It is mainly used to absorbing blood and exudates during operation and cleaning the wounds. General specification: Dia.: 3cm, 4cm, 5cm; mesh: 19×15, 26×18 EO Sterile / Non-Sterile



SH2008 Zig-zag Cotton

Absorbent cotton is pure 100% cotton fabric, degreased and bleached by advanced way. The cotton is perfect for beauty; light drainage wounds, the alcohol sterile, cleaning the medical instruments. Packing: sealed poly bag, drawstring bag or self-sealing bag General specification: 25g, 50g, 60g, 100g, 250g, 500g EO Sterile / Non-Sterile



SH2009 Cotton Ball

Soft & Pure balls are also highly absorbent, offering outstanding cleansing properties as they absorb fluids extremely quickly, making them ideal for use with toner and nail varnish or make up remover, without drips or mess.
Made from 100% pure cotton
Gentle to the skin
General specification: 0.2g, 0.3g, 0.5g, 1.0g
EO Sterile / Non-Sterile



SH2010 Cotton Applicator (Cotton Swab)

Cotton Tip Applicators are available in sterile and non sterile, different lengths and packet sizes, with a choice of wooden and plastic stems; sterile cotton tipped applicator or non-sterile cotton tipped applicators. Various size can according to the specifications of an individual purchaser, use for any skin care, peeling, facial, nail, eyebrows applications. Size: 3", 4", 6", 8" with wooden, plastic, using 100% absorbent cotton
EO Sterile / Non-Sterile



SH2401 Dental Cotton Roll

The dental cotton roll is made of 100% cotton, it is widely used on dentistry and dental cleaning.
Size: 0.4g, 0.5g, 0.7g
Size: 8 × 38mm, 10 × 38mm, 12 × 38mm
EO Sterile / Non-Sterile



SH2406 Oral Cleaning Swab

Material: PP and sponge
Color: green, red, blue, white or others;
Product Specification: 10cm, 15cm, 18cm, or as your request.
Gently cleans, moistens and refreshes mouth, teeth, gums and disinfect the wounds.
Raised foam ridges remove particles and mucous for enhanced oral hygiene



SH2409 Sponge Toothbrush

Pole material: plastic rod, length optional
Sponge colors: blue, green, pink, white, optional
Sponge shape: toothed, customizable
Sponge density: 18-30kg/m3, customizable
Packaging: Single packaging, 2 packaging, 5 packaging, 10 packaging, 20 packaging, customizable



SH2402 Saliva Ejector

Material: PVC
Length: 15cm
Packing: 100pcs/bag, 20bags/ctn, 35×29×24cm



SH2404 Dental Bib

Item: Disposable Bibs
Material: 2×20g paper Laminated 20g PE film
Color: White/green/pink/blue etc.
Packing: 125pcs/bag, 4bags/ctn 35×25×25cm



SH3213 Disposable Plastic Forceps

Material: plastic
Colour: red, blue, green, yellow, white or as customers' request.
Description: Our first and tweezers are widely used in hospitals, clinics and also in first aid.
Dimension: 10.8cm , 11.5cm, 12.5cm , 13 cm , 14cm, 18 cm etc.
Note: Logo and print can be finished according to customers' request.



SH3214 Haemostatic Forceps

Material: plastic/ABS
Colour: red, blue, green, yellow, white or as customers' request.
Dimension: 14×7.5cm



SH2101 Calico Triangular Bandages

Calico Triangular Bandages
Made of 100% cotton
Unbleached or Bleached color
With hemmed or without, with safety pins or without
Normal size: 90x90x127cm, 96x96x136cm, 100x100x140cm
Individual packing in bag



SH2102
Spandex Elastic Bandage

Material: cotton and spandex
Color: natural, flesh or bleach
Weight: 60g/m², 70g/m², 75g/m², etc.
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



SH2103
PBT Elastic Bandage

Material: cotton and polyester
Color: natural or bleach
Weight: 25g/m², 28g/m², 30g/m²
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



SH2104
Cotton Elastic Bandage

Material: 100% cotton
Color: natural or bleach
Weight: 60g/m², 70g/m², 75g/m², etc.
Size: 5cm×4.5m, 7.5cm×4.5m, 10cm×4.5m,
15cm×4.5m



SH2105
High Elastic Bandage

Material: rubber and polyester
Color: white, flesh or bleach
Weight: 80g/m², 100g/m², etc.
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



SH2106
Elastic Net Bandage

Material: rubber and polyimide
Color: white
Weight: 120g/m², etc.
Size: 0.8cm×25m, 1.7cm×25m, 2cm×25m,
2.2cm×25m, 2.6cm×25m,
3.8cm×25m, 4.4cm×25m, 5cm×25m,
6.8cm×25m, 8cm×25m



SH2108
Tubular Bandage

Material: rubber and cotton
Color: natural
Size: 4.5cm×10m, 6.25cm×10m,
5cm×10m, 6.75cm×10m, 7.5cm×10m, 8.75cm×10m,
10cm×10m, 12cm×10m, 17.5cm×10m,
21.5cm×10m, 32.5cm×10m



SH2107
Stockinette Bandage

Material: 100% Cotton
Color: Natural white, or other color color customer's demand
Introduction: Made from 100% natural cotton, which is used
for the wound and sport. It can be made into roll
Size: 5.0cm×20m, 6.0cm×20m, 8.0cm×20m,
10cm×20m, 15cm×20m



SH2109
Non-woven Self-adhesive Bandage

Material: non-woven
Weight: 35g/m², 40g/m², 45g/m², etc.
Color: 12 colors
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



SH2110
Cotton Self-adhesive Bandage

Material: cotton
Weight: 75g/m², 80g/m², 85g/m², etc.
Color: 12 colors
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



SH2112
Plaster Of Paris Bandage

Material: cotton and plaster
Weight: 390g/m², 410g/m², etc.
Color: white
Size: 5cm×2.7m, 7.5cm×2.7m,
10cm×2.7m, 15cm×2.7m



SH2115
Cast Padding

Material: natural cotton, polyester, etc
Weight: 75g/m², 80g/m², etc.
Color: white
Size: 5cm×2.7m, 7.5cm×2.7m,
10cm×2.7m, 15cm×2.7m




SH2114
Fiberglass Casting Bandage

Description: Fiber Glass casting tape consists of a lightweight knitted polyester fabric, impregnated with a water acceive polyurethane resin.
Setting time: 3-5 minutes
Weight-bearing: 20 minutes
Excellent conformability & mouldability
Tack free
100% transparent to X-rays
Very good water vapour permeability
With elasticity, and flexibility
Blue, Green, White, Pink, Red color
Size: 5cm×3.6m, 7.5cm×3.6m, 10cm×3.6m



SH2201
Zinc Oxide Plaster

Made by using cotton cloth as a base material and coating adhesive cement mainly made of zinc oxide and rubber or medical hot-melt adhesive.
Packing: simple package, plastic can, galvanized can
The product is high in stickiness, and can be cut into a saw-tooth shape
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



SH2116
Military Emergency Bandage

Size: 10cmx4m, 15cmx4m, etc
Cotton Pad Size: 10cmx10cm, 10cmx12cm, 15cmx18cm
EO Sterilization, Vacuum Package



SH2113
Cold Bandage

Cool Bandage is a cold bandage which uses a revolutionary new technique for treating injuries, muscle pain and strains. This bandage reduces recovery time from injuries by combining both the cooling and compressing process in one easy to use step.
Color: skin, blue
Size: 7.5cm×3.2m, 10cm×3.2m



SH2111
Spandex Self-adhesive Bandage

Material: spandex
Weight: 75g/m², 80g/m², 85g/m², etc.
Color: skin, or other color
Size: 5cm×4.5m, 7.5cm×4.5m,
10cm×4.5m, 15cm×4.5m



SH2202
Non-woven Surgical Tape

The nonwoven surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft, light, thin and air-permeated nonwoven cloth serving as a base material.
Packing: simple package, plastic can, with cutter
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



SH2203
PE Surgical Tape

The PE surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft PE films serving as a base material.
Packing: simple package, plastic can, with cutter
Color: transparent
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



SH2204
Silk Surgical Tape

The silk surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, the back lining material is mainly made of soft, light and thin silk cloth serving as a base material, and the adhesive is coated in a shape of mesh.
Packing: simple package, plastic can, with cutter
Color: white
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



SH2205
Sport Tape

The sports tape uses cotton cloth serving as a base material and adhesive cement mainly made of zinc oxide and rubber, and the adhesive cement is coated in a shape of strip.
Different colors available
Width: 2.5cm, 3.8cm, 5cm, etc.
Length: 10yard, 15yard, etc.



SH2209
Drilled Plaster

The product is non-toxic, low in sensitization, not irritant to the skin and subjected to EO sterilization; the product is a kind of externally applied dressing and has the effects of stopping bleeding and protecting wounds, and it is applied to affected parts such as small wounds and abrasions; and the product can prevent the intrusion of bacteria and foreign bodies and accelerate the wound healing.
Specification: PVC (natural/colored), PU, PE, elastic cloth, non-woven cloth, cartoon, etc.
Size: $\varnothing 22\text{mm}$, $\varnothing 25\text{mm}$, $38\text{cm} \times 38\text{cm}$, $72\text{cm} \times 19\text{cm}$, $76\text{cm} \times 19\text{cm}$, $76\text{cm} \times 38\text{cm}$, etc.



SH2210
Healing Capsicum Plaster

The product is mainly formed by two parts: a base material (cotton material) and separate paper. The product has the features of air permeation, moisture permeation, tight fixing, high adaptability, convenient use, etc



SH2206
Kinesio Tape

Fabric : Waterproof Elastic Fabric
Package : Simple/Individual box/Shrink/Blister
ongitudinal extension rate: 1:1.5 1:1.8
Backing paper: Silicone paper/Glassine Paper
Coating : "S"-shape adhesive coating
Color: White / Skin / yellow / Light Blue/Dark Blue / Green / Black/ Red/ Pink /Orange/Camouflage



SH2207
Wound Closure Strip

Materil: Special non-woven fabric (with a line, not with line)
Sterilization method : EO
Size : $3 \times 75\text{mm}$ / $6 \times 38\text{mm}$ / $6 \times 75\text{mm}$ / $6 \times 100\text{mm}$ / $12 \times 50\text{mm}$ / $12 \times 100\text{mm}$ / $12 \times 125\text{mm}$ / $25 \times 100\text{mm}$




SH2208
First Aid Plaster

Specification: PU, PE, Non-woven Fabric, Elastic Woven Fabric, etc.
Size: $\varnothing 22\text{mm}$, $\varnothing 25\text{mm}$, $38\text{mm} \times 38\text{mm}$, $40\text{mm} \times 10\text{mm}$, $72\text{mm} \times 19\text{mm}$, $76\text{mm} \times 19\text{mm}$, $72\text{mm} \times 25\text{mm}$, $76\text{mm} \times 38\text{cm}$, $76\text{mm} \times 50\text{mm}$, $120\text{mm} \times 20\text{mm}$, etc.



SH2212
Eye Pad

There are two kinds of Eye pad. One is non adhesive 100% Cotton eye pad. The other is made of an elastic non-woven dressing with absorbent pad. Both dressing are Specially indicated for strabismal, weak sighted therapy and nursing after ocular operation. There are different sizes for adults and children available.



SH2211

Adhesive Pad

The product is low in sensitization, excellent in biological compatibility, bacteria-free and high in air permeation, and doesn't irritate or damage the skin nor affect the natural functions of the skin.
Specification: PU, spun-laced non-woven cloth, etc.
Size: 8cm×6cm, 8cm×10cm, 10cm×10cm, 10cm×15cm, 10cm×20cm, 10cm×25cm, etc.



Sh2306

Vaseline Gauze Swab (Paraffin Gauze)

The product is made from medical degreased gauze or non-woven together with Vaseline. It can lubricate the skin and protect the skin from cracks. It is widely used on clinic.



SH2311

Sting Relief Medicated Pad

Contents: Antipruritic solution
Non-Woven Size: 65×30mm, 65×58mm, etc.
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm



SH2307 Wet Wipes

Specifications:

Material: Spunlace Nonwoven, ES hot air through Nonwoven, Thermal-bonded Nonwoven, Air-laid paper, Wet Strength Paper, Bamboo Fiber.
Weight: 25-80gsm for your choice
Size: 15x20cm, 16x20cm, 18x20cm, 20 x 30cm, etc.
Wipes count: 1 to 200pcs for your choice.
Folding Size: Z fold, M fold.
Formula: Gentle and safe, Alcohol Free, all kinds of fragrance for your choice.
Fragrance: Lemon, Lavender, Camellia, etc.
Function extracts for your choice, such as Aloe Vera & Vitamin E and so on.
Packaging type: PE/PET film, plastic canister/tube, or plastic box, etc.



SH2301

Alcohol Prep Pad (Alcohol Swab)

Contents: 70% Isopropyl Alcohol or 70% Isopropyl Alcohol + Chlorhexidine
Non-Woven Size : 65×30mm, 65×58mm, etc.
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm



SH2302

Povidone-Iodine Prep Pad

Contents: 10 % Povidone-Iodine Solution
Non-Woven Size : 65×30mm, 65×58mm, etc.
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm



SH2308

Alcohol Wipes

Contents: 70% Isopropyl Alcohol or 75% Ethyl Alcohol
Size : 15x20cm, 20x20cm, 20x30cm, etc.
Material: Spunlace nonwoven, ES hot air through nonwoven, thermal-bonded nonwoven, air-laid paper, wet strength paper, bamboo fiber
Wipes Count: 1pcs, 10pcs, 20pcs, 50pcs, 100pcs

SH2309 Disinfectant Wipes

Specifications:

Material: Spunlace Nonwoven, ES hot air through Nonwoven, Thermal-bonded Nonwoven, Air-laid paper, Wet Strength Paper, Bamboo Fiber.
Weight: 25-80gsm for your choice
Size: 15x20cm, 16x20cm, 18x20cm, 20 x 30cm, etc.
Wipes count: 1 to 200pcs for your choice.
Folding Size: Z fold, M fold.
Fragrance: Lemon, Lavender, Camellia, etc.
Packaging type: PE/PET film, plastic canister/tube, or plastic box, etc.
Different agent solution (680g, 800g, etc.) for choice each loop agent solution contain 25g ethanol (94%) and 35g propan-1-ol.



SH2304

Povidone-Iodine Swabstick

Contents: 10 % Povidone-Iodine Solution
Swabstick Size: 70mm, 100mm
Packet Size : 110×40mm, 140×40 mm



SH2305

Alcohol Swabstick

Contents: 70% Isopropyl Alcohol or 70% Isopropyl Alcohol + 2% Chlorhexidine
Swabstick: 70mm, 100mm
Pack Size: 110×40mm, 140×40mm

MEDICAL NON-WOVEN & ACCESSORIES



SH2601
Round Nurse Cap

Material: non-woven PP
Weight: 10g/m2, 12g/m2, 14g/m2, 16g/m2, etc
Size: 18", 21", etc
Color: white, blue, green, etc



SH2602
Clip Cap

Material: non-woven PP
Weight: 10g/m2, 12g/m2, 14g/m2, etc
Size: 20", 21", 22", etc
Type: single elastic/double elastic
Color: white, blue, green, etc



SH2603
Doctor Cap

Material: non-woven PP
Weight: 20g/m2, 25g/m2, 30g/m2, etc
Size: 61×14cm, 64×23cm, etc
Type: tie/elastic back
Color: white, blue, green, etc



SH2605
Surgical Hood

Material: non-woven PP
Weight: 16g/m2, 20g/m2, 25g/m2, 30g/m2, etc
Size: 30×124cm, etc
Type: tie
Color: white, blue, green, etc



SH2610
Non-woven Face Mask

Material: non-woven PP or non-woven PP and filter paper
Specification: ear-loop/tie and 2ply/3ply/4ply
Size: 17.5×9.5cm, 14.5×9cm
Color: white, blue, green, etc



SH2611
Non-woven Face Mask with Charcoal
Description:
Adjustable nose seal,easy to fold and carry. Mini-fiber high electrostatic filter, filtration efficiency of dust & powder is over 95%. Repel continually contamination of dangerous blood and saliva.
Specification: Ear-loop type: 4 ply
Detailed weight: 4 ply: 25g/m²+20g/m²+40g/m²+25g/m²
Size: 18×9cm
Color: white, green, blue etc.
Packing: 50pcs/box, 40boxes/carton, 520×380×300cm



SH2612 Paper Face Mask
Material: high quality wet-strength paper
Specification: ear-loop and 1ply/2ply
Size: 20 × 7cm
Color: white



SH2614
Face Mask With Shield

Composition: non-woven face mask and PET shield (prevent mist)
Specification: ear-loop and 3ply/4ply
Size: 17.5x9.5cm
Color: white,blue,green,etc



SH2615
Face Guard

Material: PET
Type: elastic double-sided anti-fog
Size: 320x220mm
Thickness: 0.25mm
Packing: 1pc/pe bag



SH2618
Protective Spectacle

Impact resistant PC lens, Lens thickness: 2mm
Integrated side legs
High strength structure
Integrated wide side shields



SH2619
Protective Goggles

Material: PVC+PC
Clear, Anti-fog
Packing: 1pc/poly bag



SH2613
Simple-Dust-Proof FaceMask

Conforms to FFP1 & FFP2
Fiber glass Free, Adjustable contoured aluminum nosepiece helps to ensure proper fit.
Non-woven PP material: 16g, 18g, 20g, etc.
Color: white, blue, etc.



SH2620
New-Type-Dust-Proof Face Mask

Material: terylene (conforms to American NEL SON N95 test standards)
Specification: ear-loop
Size: 220g,240g,260g,etc
Color: white,blue,etc



Sh2616
N95 Face Mask (with/without valve)

Material: polypropylene (conforms to American NEL SON N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white



SH2617
Folded Facemask

Material: polypropylene (conforms to American NEL SON N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white

SH2621

Surgical Gown

Standard and Reinforced
Material: Spunlace, SMS, SU/US nonwoven and so on
Gram weight of material: from 35g/m2 to 50g/m2
Color can be decided by customer
Whether sterilize is be decided by customer
Packaging way can be bulk or individual
Size: S, M, L, XL, XXL, XXXL



SH2622

Disposable Isolation Gown

Material: SPP or PP+PE non-woven fabrics
or as customers' requirement
Color: Blue, yellow, green and so on.
Size: S/M/L/XL/XXL
Details: knitted or elastic cuff
with velcro or ties on the neck and ties on the waist.



SH2623

CPE Coat

Material: CPE
Color: Blue
Size: 62g, 125×200cm, etc
Characteristics:
prevent and isolate dust,
particle, alcohol, blood
Bacteria and virus invading



SH2621R

Reinforced Surgical Gown

Specification:
Material: Spunbond Meltdown Spunbond (SMS)
Color: white/yellow/green
Weight: 40gsm/45gsm/48gsm/50gsm/55gsm etc.
Size: S/M/L/XL/XXL
Suited to hospital use & for surgery
With ties at back of the neck and waist
Comfortable sleeve tailoring for maximum
Freedom of movement
Soft knit cuffs for wrist comfort
Secure sleeve seams for protection
Extra protection reinforce on sleeve & chest
And full-body on selected gowns
Spunbond meltdown spunbond(SMS) material gown breathable,
quiet and supple,
High repellency, AAMI standard, sterile with comfortable design
for operating theater
With knitted cuff, Velcro neck back closure, side waist band
Reinforced arm and additional liner at front
1 hand towel is available
1 wrap towel is available
Sterile method: non-sterile or EO sterile
Packing: 10pcs/poly bag or 1 pc/sterile bag



SH2624

Lab Coat

Material: non-woven PP/PE
Color: white/blue/yellow/green
Weight: 30g, 35g, 40g, 45g, etc
Size: M/L/XL/XXL



SH2626

Disposable Protective Coverall

Chemical protective coverall includes: liquid-tight
protective clothing , spray-tight protective clothing
protectiveclothing for use against solid particulates
protective clothing against liquid chemicals and so on.
Material: PP/PP+PE/SMS and so on
Weight: 35-60g
Size: M/L/XL/XXL
Color: Blue/white/yellow/ dark blue
Application: provide reliable protective products for
industry/lab/ nuclear radiation/ agriculture/ military
combat readiness and so on



SH2625

Disposable Scrub Suit

Size: S/M/L/XL/XXL and so on
Material: SS/PP/SMS 20g - 60g
Color: Blue/white/yellow/ dark blue
Packing: 1pc/PE bag, 50bags/carton



SH2632 PE Shoe Cover

Material: PE, with comfortable elastic opening
Color: blue, green, etc
Weight: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
Size: 36×15cm, 40×16cm
Thickness: 0.018mm, 0.02mm, 0.021mm and so on



SH2633 CPE Shoe Cover

Material: CPE, with comfortable elastic opening
Color: blue, green, etc
Weight: 1.6g, 1.8g, 2.0g, 2.1g, 2.5g
Size: 36×15cm, 40×16cm
Thickness: 0.018mm, 0.02mm, 0.021mm and so on



SH2635 Non-woven Shoe Cover

Material: non-woven PP
Color: white, blue
Regular/Non-skid
Weight: 30g/m², 40g/m², 45g, etc
Size: 36×15cm, 38×18cm, 40×16cm



SH2631 PE Sleeve

Made of PE, with comfortable elastic
Color: white, blue and red, etc
Thickness: 0.018mm
Size: 20×40cm, 22×42cm, 22×46cm



SH2636 PE Apron


Made of polyethylene, with adjustable tie
Color: white, blue, etc
Size: 110×70cm, 135×71cm, 160×90cm
Thickness: 0.018mm, 0.02mm and so on



SH2641

Non-woven Bed Sheet

Soft touch and harmless to your skin
100%biodegradable and environmentally friendly
SPEC: light blue , 60×180cm
30gsm other size and specification are also acceptable



SH2642

Non-woven Bed Cover

Disposable bed cover made of pp non-woven fabric
helps you to be clean & hygiene on your travel and
regular purpose
SPEC: 30g, 210×90×20cm, other size and specification
are also acceptable



SH2643

Non-woven Pillow cover

Disposable pillow covers are used widely in hotels
hospitals, home, and airline and in all travels


Specification:
80gsm,50×60cm,with 5cm gusset,
other size and specification are also acceptable
Packing: 100pcs/polybag, 1bag/ctn
Ctn size: 54×32×38 /100pcs



SH2651

Under Pad

Material: non-woven fabrics+2ply medical paper
+ fluff pulp+PE film, non-woven fabrics+4ply
medical paper+PE film, etc
with or without SAP
Size: 40×60cm, 60×90cm, etc
Weight can be decided by customer



SH2652

Adult Diaper

First layer: Non-woven fiber, natural soft
and breathable.
Second layer: fuff pulp+super absorbent polymer
+classic tissue.
Breathable PE film with cloth like back Non-woven
or casting film.
Standard Size:M: 650×800mm
L: 750×900mm, XL:840×1030mm




SH3211 Surgical Brush

Medical Brush has an effective, soft and skin-friendly brush
head for the nails, knuckles and skin.
Attached to the brush is a sponge designed for cleaning the
hands and forearms.
In addition a separate nail cleaner is supplied in the packaging.
The Medical Brush saves on handling costs since brush, soap
and nail cleaner are all supplied in one packaging.
In addition to a dry sterile version, one has a choice of several
impregnated versions:
Medical Brush Dry Sterile
Medical Brush Chlorhexidine
Medical Brush PVP Iodine



SH3212 Surgical Hand Brush

Surgical scrub brush is used to significantly reduce the number
of micro-organisms
on the hands and forearms before performing surgeries.
The nail cleaner feature makes removal of micro-organisms
under the nail quick and easy
All kinds of spray cans are available



SH2701 Universal Pack

- 1 x Back Table Cover 140x190cm
- 1 x Mayo Stand Cover 78x145cm
- 2 x Adhesive Side Drape 75x90cm
- 1 x Adhesive Bottom Drape 150x180cm
- 1 x Adhesive Top Drape 150x240cm
- 4 x Hand Towel 30x40cm
- 1 x OP Tape Strip 10x50cm
- 2 x Reinforced Surgical Gown, Lsize



SH2702 Disposable Sterile C-Section Surgical

- 1 x Back Table Cover 140x190cm
- 1 x Mayo Stand Cover 76x145cm
- 1 x C-section Drape, fluid collection pouch with 2 suction port 72x80cm and forming wire, full incision film 40x40cm, tube holders, Armboard covers, 240x300x195cm
- 4 x Hand Towel 30x40cm
- 3 x Reinforced Surgical Gown with 2pcs hand towel, XL
- 2 x Baby Blanket 60x90cm



SH2703 Disposable Sterile Angiographics

- 1 x Trolley Cover 140x150cm
- 1 x Femoral Angiography Drape 240x300cm, femoral Fen. Φ 12cm, incision film
- 1 x Equipment Protector Cover 80x110cm
- 1 x Image Intensifier Cover, 140cm20x Gauze, 8ply 10x10cm .
- 1 x Sponge
- 1 x Scalpel #11
- 2 x Surgical gown XL
- 1 x Plastic Bowl 500ml
- 1 x Plastic Cup 60ml
- 1 x Tray 30x30x5cm



SH2704 Knee Arthroscopy Pack with Gown

- 1x Back Table Cover 140x190cm
- 1x Mayo Stand Cover 78x145cm
- 1x Arthroscopy T Drape 220x280x290cm, Fluid collection pouch 83x78cm, Circular Elastic fenestration 6cm, Absorbent reinforcement 75x120cm
- 1x Suture Bag
- 2x OP Tape Strip, 10x50cm
- 4x Hand Towel 30x40cm
- 2x Reinforced Surgical Gown with 2pcs hand towel, XL



OTHERS



SH2503 Hot Water Bag

Material: calcium carbonate \geq 30% rubber
Size: 500cc, 1000cc, 1750cc, 2000cc, etc
Packing: 1pc/PE bag




SH2512 Ear Syringe

Material: calcium carbonate \geq 25% rubber or PVC
Size: 30ml, 60ml, 90ml, etc
Packing: 1pc/PE bag



SH250602 Hot/Cold Pack

Material: PVC+carbamide + distilled water
Weight: 40g
Size: 9.8 \times 8cm, 11.5 \times 11cm, etc.



SH250603 Hot/Cold Pack

Material: PVC+carbamide + distilled water
Weight: 65g
Size: 9 \times 9cm, 10 \times 10cm, 12 \times 12cm, etc.



SH2511 Air Cushion

Material: calcium carbonate \geq 30% rubber
Size: ϕ 35cm, ϕ 40cm, ϕ 42.5cm, ϕ 45cm, etc
Packing: 1pc/PE bag



SH2504 Fabric Ice Bag

Made of canvas
Size: 6", 9", 11", etc
Packing: 1pc/PE bag or 1pc/paper box



SH3301 Cupping Massage

Gift Package



SH2902 Umbilical Cord Clamp

Material: ABS, PE, PP
Color: blue, white, etc
Size: 5cm, 5.5cm



SH2505 Instant Ice Bag

Material: PE+PA infilling (carbamide + distilled water)
Weight: 160g, 280g
Size: 19.5 \times 12cm, 24 \times 12cm



SH250601 Hot/Cold Pack

Material: nylon taffeta infilling (glycerol)
Weight: 200g
Size: 20.5 \times 15.5cm



**SH2903
Umbilical Cord Clamp Cutter**



**SH2904
Amniotomy Hook**




SH3101 Heat-Sealing Flat Pouch

Size: 57×130mm , 70×230mm
90×230mm , 134×255mm
190×330mm



SH3102 Self-Sealing Flat Pouch

Size: 50×130mm , 90×230mm
135×255mm , 190×330mm
250×370mm



SH3108 Medical Crepe Paper

Size: 30×30cm
50×50cm
70×70cm
120×120cm etc.



**SH3105
Sterilization Indicator Tape**

Model: EO, Autoclave, Plasma and so on



SH3104 Paper/Paper Gusseted Pouch

Size: 90×50×150mm , 90×50×250mm
110×30×190mm , 140×75×250mm etc.




SH3103 Heat-Sealing Gusseted Pouch

Size: 50mm×250mm , 150mm×200mm
100mm×250mm , 140mm×250mm
150mm×300mm , 200mm×350mm



**SH3106
Steam Indicator Strip**

Size: 15×100mm, 15×200mm, 200pcs/box, etc.



**SH3003
Disposable E.C.G Electrode**

Model: non-woven, foam, micropore, medical tape, etc



SH3101R Heat-Sealing Flat Reel

Size: 55mm×200m , 75mm×200m
100mm×200m , 150mm×200m
200mm×200m




SH3103R Heat-Sealing Gusseted Reel

Size: 50mm×100m , 55mm×100m
75mm×100m , 100mm×100m
150mm×100m



**SH3004
Electrosurgical Pad**

Size: Child, Adult
Material: PE Foam or Non-woven, Aluminum-foil, Hydro-gel
Packing: 1pc/pouch



SH3005 E.C.G Paper

Classification : single lead electrocardiograph
three lead electrocardiograph
six lead electrocardiograph
twelve lead electrocardiograph



SH3007 Disposable Oxi-Sensor

7 pin/9 pin, Neonate/Adult
1pc/blister packing, 40pcs/box. 4boxes/ctn



SH0651 Tourniquet

Highly flexible,
latex free cotton ribbon which is gentle on the skin.
sturdy snap-on clip
quick and safe releasing



Sh06521 Pre-cut Tourniquet

Material: TPE
Size: 1"x18" , 0.6mm in thick, with/without texture
1pc/roll or Flat-Packing



SH06522 Pre-cut Roll Tourniquet

Material: TPE
Size: 1"x18" , 0.6mm in thick
25 each per big roll , per box



SH2907 I.D. Bracelet

Material: PVC
Model: write-on, insert-card
Size: infant, adult, adult/infant
Color: blue, pink, white and so on



SH2908 Eyecare Phototherapy Mask

Size: N-A Head Circumference 30-38cm
N-B Head Circumference 24-33cm
N-C Head Circumference 20-28cm
T Head Circumference 24-35cm

Sharps Container



SH280101

Capacity: 0.8L
Product Size:
Top Dia: 125mm
Bottom Dia: 95mm
Height: 115mm



SH280102

Capacity: 3.0L
Product Size:
Top Dia: 170mm
Bottom Dia: 150mm
Height: 150mm

SH280103

Capacity: 3.0L
Product Size:
Top Dia: 200mm
Bottom Dia: 160mm
Height: 140mm



SH280104

Capacity: 3.0L
Product Size:
Top Dia: 200mm
Bottom Dia: 120mm
Height: 170mm



SH280105

Capacity: 1.0L
Product Size:
Top Dia: 110mm
Bottom Dia: 90mm
Height: 160mm





SH280108

Capacity: 3.0L
Product Size:
Top Dia: 150mm
Bottom Dia: 120mm
Height: 260mm



**SH280112
Needle Cutter**

Item: Needle Cutter
Material: ABS + Blade
Description:
Dia 11.33cm x 9.6cm, manual



**SH3201
Dressing Change Kit**

Product Composition :
1.povidine-iodine cotton ball
2.dry cotton balls
3.plastic tweezers
4.PE glove
5.drape
6.ahesive wound dressing
7.tray



**SH3202
Suture Remove Kit**

Product Composition :
1.drape
2.PE glove
3.povidine-iodine cotton ball
4.suture scissor
5.plastic tweezers
6.tweezers (metal)
7.tray
8.gauzes



**SH280109
Paper Safety Box**

Capacity: 5L, 10L etc.



**SH2802
Biohazard Specimen Bag**

Size: 8"x10"
Thickness : 2mil 3ply
Size: 6"x10"
Thickness: 2mil 3ply



**SH3203
Oral Care Kit**

Product Composition :
1.bib
2.probe
3.tweezers (metal)
4.tray
5.moatn mirror
6.dry cotton balls



**SH3203
PICC Puncture Care Kit**

Product Composition :
1.wrapping paper
2.drape
3.transparent dressing
4.latex gloves
5.povidine-iodine prep swab
6.alcohol swab
7.wound closure
8.non-woven gauze
9.alcohol preppad