



SeaSight®

NINGBO SEASIGHT MEDICAL INSTRUMENTS CO., LTD.

Office Add: 9F, Ningshing Building, No.666 Jiangnan Road, Ningbo, China.
Factory Add: No.508, Yindong Road(N), Yinzhou Economic
Development Zone, Ningbo315145, China
Tel: 0086-574-81877107 Fax: 0086-574-81877057
<http://www.seasight-med.com>



SEASIGHT
Medical
Instruments
always by your side

ABOUT US

»»»» OUR COMPANY



SeaSight®

always by your side

Establishment

SEASIGHT is established in the year 2003 as a professional manufacturer and supplier for medical and lab products in China. Our 1st factory- MFLAB Medical with production area 15,500 m2. It includes 100000 grade purification workshops and cleaning rooms, totally 9,200 m2. Our 2nd Factory-SEASIGHT Medical (production area covering 52,000 m2) will open soon.

Team Work

SEASIGHT has a high efficient and professional team, full of wisdom and passion, with more than 20 years of working experience. We take the diligent, professional and responsible spirit, committed to provide customers good products and service.

Export Reputation

Our SEASIGHT products have been exported to over 150 countries worldwide through a network of distributors and subsidiaries, covering North American, South American, European, Asian and African countries.

Leading Products

SEASIGHT has a wide range of products in business with leading products respiratory products, disposable medical products, medical equipment and laboratory products. We can supply clients with OEM or ODM or OBM.

Brands

There are Seven brands in our group: Seasight, Sunflos, Medland Care, Foyocare, Foyomed, FoyoLab and Mflab in distribution.

Quality System

Our factory MFLAB Medical audited by TUV Rheinland every year, more than 40 products has got CE certificates issued by TUV Rheinland, and also passed by GMP, ISO13485, US FDA on-site Audit, Saudi Arabian FDA on-site Audit, Malaysian MOH on-site Audit, and Korea MOH on-site Audit. We have also got a number of patents.

Principle

We SEASIGHT adheres to "Always by your side" principle to strive to achieve excellent quality, competitive price and first-rate service in operation, and cherish every business opportunity.



LB1310

Oxygen Mask With Reservoir Bag

Made of clear, medical grade PVC.
With reservoir bag.
With adjustable nose clip, elastic strap with/without latex.
Oxygen mask with or without check valve.
Color: green transparent, white transparent and light blue transparent are available.
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Cat.No.	Size
LB131000	XS(Neonate)
LB131001	S(Pediatric Standard)
LB131002	M(Pediatric Elongated)
LB131003	L(Adult Asian)
LB131004	L+(Adult European)
LB131005	XL(Adult Elongated)




LB2011

Nebulizer Mask

Made of clear, medical grade PVC.
Operates even at horizontal position capacity: 6cc/10cc.
With adjustable nose clip, elastic strap with/without latex.
Color: green transparent, white transparent and light blue transparent are available.
Length: with 2m tubing, different length can be customized,
Anti-crush tube can ensure oxygen follow even if the tube is kinked.
Packing: packed in individual PE pack. 50 pcs/carton or 100 pcs/carton.

Cat.No.	Size
LB201100	XS(Neonate)
LB201101	S(Pediatric Standard)
LB201102	M(Pediatric Elongated)
LB201103	L(Adult Asian)
LB201104	L+(Adult European)
LB201105	XL(Adult Elongated)




LB6100
Suction Yankauer Handle

LB6201



LB6203



Size	Type
LB6101	Crown tip with vent
LB6102	Crown tip without vent
LB6103	Flat tip with vent
LB6104	Flat tip without vent

LB6202



LB6200
Connector

Supplied with universal connector,
can connect with smooth suction yankauer handle.

Size	Type	Color
LB6201	Long Connector	Blue, Green or Transparent
LB6202	Short Connector	Light Blue or Navy Blue
LB6203	Trumpet Connector	Transparent or Blue

Connecting Tube

Made of transparent material for better visualization
Striated walls of tube provide superior strength
and anti- Kinking

Size	ID(mm)	OD(mm)
1/4"	6.5±0.2	9.7±0.2
3/16"	4.7±0.2	8.2±0.2




LB60
Connecting Tube without Yankauer Handle



LB61
Connecting Tube with Yankauer Handle

LB6204
Optional Accessories



Made from Non-toxic PVC, clear and soft
Large Lumen resists clogging and the transparency
Allows Clear visualization of liquids
Crown tip, with/without Vent or plain tip
with/without Vent.
Packing: Packed in individual blister
packing or Double packing, 20pcs/carton.
Length can be customized.

Circuits

Seasight offers a wide range of anesthesia circuits available in different lengths and configurations for your clinical needs and requirements, including different tubing types: corrugated, smoothbore, extendable, coaxial, duo limbo and different tubing sizes: adult 22mm, pediatric 15mm.

Corrugated Circuits

- ISO standard connectors, matched with anesthesia and respiratory machines
- Channel for medical gases (anesthesia gas oxygen) entering into patients
- Flexible, soft, transparent, airtight
- Good compliance low airflow resistance



Extendable Circuits


With adjustable length, the extendable circuit allows us to optimize the compression volume as required, which is easier for transportation storage and cleaning.



Smoothbore Circuits

Smoothbore circuit unifies a soft connector for the machine end with another hard connector to improve ventilatory support.

- Advantages:
- Unique swivel tube, preventing breath hold
 - Smooth inner wall, not accumulating water,improving safety.



Good Elasticity, Good Airtightness / For Anesthesia Room and ICU



Anesthesia Circuits



Coaxial Circuits



Duo Limbo Circuits



Breathing Circuits



Non-invasive Circuits



Circuits Accessories

always by your side

E-mail:market@seasight-med.com
www.seasight-med.com



LB1810 Tri-ball Incentive Spirometer

600cc/sec+900cc/sec+1200cc/sec

The Tri-ball Incentive Spirometer is designed to encourage the user to take deep breath. The spirometer provides patient feedback by showing how well the patient is expanding their lungs. This deep breathing encourages the clearance of phlegm from the lungs.

LB1710 Peak Flow Meter

It is used to detect the pulmonary for the asthma patients

Material: ABS

Size: Adult/Child

Packing: 1pc/box, 100pcs/carton



LB1720 Peak Flow Meter

It is used to detect the pulmonary for the asthma patients

Material: ABS

Size: 30/50 x 155mm

Packing: 1pc/box, 200pcs/carton





LB2120

Spacer For Aerosol

Cat.No.	Size	Applicable People	Color	Mask Weight
LB212001	S	0-18 Months	Red	19.7G
LB212002	M	1-5 Years	Yellow	23.7G
LB212003	L	5+ Years	Purple/blue	43.8G

- Spacer Body Size** : 50x153 (bottom diameter x height)
- Volume** : 175ml
- Mouthpiece** : 22.5x15.5mm(top of mouthpiece dimension)
- Spacer Body Weight**: 50g
- Material** : 100% Latex free, Silicone Rubber, Anti-static PETG and TPR.

LB3110

Silicone Resuscitator

LB3110L

Liquid Silicone Resuscitator

Silicone resuscitator is intended for pulmonary resuscitation
The set includes oxygen reservoir bag, silicone mask and oxygen tubing
Reusable, autoclave 134℃ helping to prevent cross infection and contamination
Guaranteed for 5 years of alternatively 20 times autoclave cycles
Patient connector is 22/15mm.

Cat.No. Silicone	Cat.No. Liquid Silicone	Size	Specification	Package
LB311001	LB311001L	Adult: Adults and Children	Reusable adult 1500ml resuscitator(W/60cm H ₂ O POP-OFF) Silicone mask 5# 2000ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:60cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available
LB311002	LB311002L	Pediatric: Children and Infants	Reusable pediatric 550ml resuscitator(W/40cm H ₂ O POP-OFF) Silicone mask 2# 1600ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:40cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available
LB311003	LB311003L	Infant: Infants and Neonates	Reusable infant 280ml resuscitator(W/40cm H ₂ O POP-OFF) Silicone mask 1# for infant Or silicone mask 0# for neonate. 1600ml oxygen reservoir bag Oxygen tubing 2.1meter Pressure relief valve:40cm H ₂ O pressure relief	1pc/pack, PE bag, Paper box, PP box available



LB3210
Disposable Silicone Laryngeal Mask

The tube is made from 100% medical silicone. No poison and harmless
Sterilization: EO
Packing: 1pc/blister packing

Cat.No.	Size	Applicable People
LB321010	1#	<5kg
LB321015	1.5#	5-10kg
LB321020	2#	10-20kg
LB321025	2.5#	20-30kg
LB321030	3#	30-50kg
LB321040	4#	50-70kg
LB321050	5#	>70kg

Standard Endotracheal Tube (Oral/Nasal Endotracheal Tube)

Endotracheal tube is a temporarily breathing device that is inserted into patient's trachea in order to maintain an open air passage and prevent aspiration of mucus and gastric contents. It is designed for general anaesthesia , intensive care and emergency case for airway management and mechanical ventilation.



LB5010C
High Volume, Low Pressure Cuff

Cat.No.	Size ID(mm)	Cat.No.	Size ID(mm)
LB501025C	2.5	LB501065C	6.5
LB501030C	3.0	LB501070C	7.0
LB501035C	3.5	LB501075C	7.5
LB501040C	4.0	LB501080C	8.0
LB501045C	4.5	LB501085C	8.5
LB501050C	5.0	LB501090C	9.0
LB501055C	5.5	LB501095C	9.5
LB501060C	6.0	LB501100C	10.0




LB5010
Uncuffed

Cat.No.	Size ID(mm)	Cat.No.	Size ID(mm)
LB501025	2.5	LB501065	6.5
LB501030	3.0	LB501070	7.0
LB501035	3.5	LB501075	7.5
LB501040	4.0	LB501080	8.0
LB501045	4.5	LB501085	8.5
LB501050	5.0	LB501090	9.0
LB501055	5.5	LB501095	9.5
LB501060	6.0	LB501100	10.0



LB6200
Suction Tube

Made of medical grade PVC.
Color-coded for size identification.
Open distal end with 2 side holes.
Different connector can be chosen.
Standard with 50cm length.
Different length is available.



Size ID (Fr)	Cat.No.:LB6210 (Plain Type)	Cat.No.:LB6220 (Cap-Cone)	Cat.No.:LB6230 (Y-Tpye)	Cat.No.:LB6240 (T-Tpye)	Cat.No.:LB6250 (Y-Tpye Transparent)	Color
Fr6	LB621006	LB622006	LB623006	LB624006	LB625006	Light Green
Fr8	LB621008	LB622008	LB623008	LB624008	LB625008	Blue
Fr10	LB621010	LB622010	LB623010	LB624010	LB625010	Black
Fr12	LB621012	LB622012	LB623012	LB624012	LB625012	White
Fr14	LB621014	LB622014	LB623014	LB624014	LB625014	Green
Fr16	LB621016	LB622016	LB623016	LB624016	LB625016	Orange
Fr18	LB621018	LB622018	LB623018	LB624018	LB625018	Red
Fr20	LB621020	LB622020	LB623020	LB624020	LB625020	Yellow
Fr22	LB621022	LB622022	LB623022	LB624022	LB625022	Purple

Different Connector is available.

Latex Foley Catheter

Intended Use /Indications

The Foley Catheter is passed through the urethra during urinary catheterization and into the bladder to drain urine, or for inserting fluids into the bladder. It is used in departments of urology, internal medicine, surgery, obstetrics, and gynecology for drainage of urine and medication. It is also used for patients suffering from moving with difficulty or being completely bed-ridden.



Size (Fr/Ch)	Balloon (cc)	Color Code	Length (mm)	Cat.No.
2-way pediatric				
6	3	Light red	270	SH010106
8	5	Black	270	SH010108
10	5	Grey	270	SH010110
2-way standard				
12	30	White	400	SH010112
14	30	Green	400	SH010114
16	30	Orange	400	SH010116
18	30	Red	400	SH010118
20	30	Yellow	400	SH010120
22	30	Violet	400	SH010122
24	30	Blue	400	SH010124
26	30	Pink	400	SH010126
2-way female				
12	30	White	270	SH010312
14	30	Green	270	SH010314
16	30	Orange	270	SH010316
18	30	Red	270	SH010318
20	30	Yellow	270	SH010320
22	30	Violet	270	SH010322
3-way standard				
16	30	Orange	400	SH010216
18	30	Red	400	SH010218
20	30	Yellow	400	SH010220
22	30	Violet	400	SH010222
24	30	Blue	400	SH010224
26	30	Pink	400	SH010226

Features


- ★ Made from natural latex
- ★ Silicone coated
- ★ Color coding for size identification
- ★ Sterilized by EO
- ★ For single use only
- ★ Available with different types
- ★ Available with rubber&plastic valve

Natural Latex Rubber Condoms



Specifications
Normal Width: 49 + / - 2mm, 52 + / - 2mm, 55 + / - 2mm
Standard: ISO4074-2015
100% pure natural latex concentrated and compounded
100% electrically tested to ensure reliability CE & FDA approved

Material: Natural Latex
Odour: Banana, Orange, Cherry, Strawberry, Mint, Pear, Apple, Grape
Different Styles: Classic Plain Condoms
Ribbed Condoms
Super Thin Condoms
Dotted Condoms
Extra Strong Condoms
3 In 1 Condoms
Super Dotted Condoms





SH0802 Surgical Gloves

Size: 6.0#, 6.5#, 7.0#, 7.5#, 8.0#, 8.5#, 9.0#
Powdered/Powder free
1pair/pouch, gamma sterile
AQL: 1.5
Packing: 50pairs/box, 500 pairs/ctn

SH0803 Latex Examination Gloves

Size: S, M, L
Powdered/Powder free
100pcs/box, non-sterile
AQL: 2.5
Packing: 100pcs/box, 1000 pcs/ctn



SH0804 Nitrile Examination Gloves

Size: S, M, L
Color: Blue, Black, Lavender. etc
Powder free
100pcs/box, non-sterile
AQL: 1.5/2.5/4.0
Packing: 100pcs/box, 1000 pcs/ctn

SH0805 Vinyl Gloves

Size: S (3.0g, 3.5g, 4.0g) ; M (3.5g, 4.0g, 4.5g) ;
L (4.0g, 4.5g, 5.0g)
Color: Blue, White
Powdered/Powder free
100pcs/box, non-sterile
AQL: 4.0
Packing: 100pcs/box, 1000 pcs/ctn

SH0810 Vaginal Speculum



SH0811 Side Screw Type

Size: S,M,L

SH0812 Fastener Type

Size: S,M,L



SH0813 Spanish Type

Size: S,M,L

SH0814 Middle Screw Type

Size: S,M,L



SH0242
Closed Wound Drainage System (spring)

Our closed wound drainage system, including 400cc, 3-spring evacuator, PVC tubing and "Y"connector, and drain with trocar. The 3-spring evacuator design is minimize collapse-resistance, packed in double liner Tyvek, to make sure medical cleanness on O.R.operation.
Capacity: SH0242 (200ml, 400ml,800ml)
Silicone or PVC tube, 120cm Length
Sterilization: EO
Packing: Individual blister Packing, 20pcs/ctn



SH0235
Closed Wound Drainage System Kit

Made of 100% imported medical grade silicone.
This product includes 100cc/200cc/400cc Silicone Reservoir, 7mm/10mm Flat drains or other round drains and stainless steel trocar.
Can use separately or using Trocar for insertion and combining with silicone Reservoir for fluids syction and collection.
Used for various surgeries.
Supplied sterile, EO Sterilization.



SH0236
Silicone Reservoir/Bulb

Made of 100% imported medical- grade silicone,with anti- reflux valve.
Can compatible with any of following drains:
Flat Drain, Fluted Drain, Round Drain and Perforated Drain.
Squeeze the reservoir and connect with drains for suction, the low-vacuum reservoir can make stable and effective suction for exudates.

Cat. NO.	Size(cc)	Marked Position(ml)
SH023601	100	25
SH023602	200	30
SH023603	400	40

*Clamp available.

SH0604
Safety Clip Syringe

1ml, 2/3ml, 5ml, 10ml, 20ml
Safety syringe with protective cap
Luer lock
After injection, push the clip over the needle and lock it onto the clip base
Sterilization: EO



SH0610
Insulin Pen Needle

Ultra fine cannula for maximum patient comfort, patient point with triple beveling for easy and pain decreases
The needle is fit for Humapen (EliLily, USA)
NovopenIII (Novo Nordisk, Denmark)
Sterile, non-pyrogenic
Sterilization: EO



SH0911
Pressure Activated Safety Lancets

Needle is fully shielded before and after use
Pre-loaded, pressure activated design
Ensures a good blood sample every time
Ultra sharp tri-bevel needle or bi-bevel blade, high speed penetration
Variety of needle specifications to suit most capillary blood requirements
Sterilization: Gamma
Size: 21G, 23G, 26G





SH2001 Gauze Swabs

The pads are folded all by machine, pure 100% cotton yarn ensure the product soft and adherent. Superior absorbency makes the pads perfect for absorbing blood and exudates.
In accordance with the customers' requirements we can produce different kinds of pads, such as folded and unfolded, with x-ray and non x-ray
The adherent pads are perfect for operation.
General specification: 2"×2", 3"×3", 4"×4"
PLY: 8ply, 12ply, 16ply
Mesh: 19×10, 19×15, 20×12, 26×18
EO Sterile / Non-Sterile



SH2202 Non-woven Surgical Tape

The nonwoven surgical tape mainly uses low-sensitization medical hot-melt adhesive, the back face is coated with an uncoiling agent, and the back lining material is mainly made of soft, light, thin and air-permeated nonwoven cloth serving as a base material.
Packing: simple package, plastic can, with cutter
Color: white, flesh
Width: 1.25cm, 2.5cm, 5cm, 7.5cm, 10cm, etc.
Length: 4.5m, 5m, 10yard, 10m, etc.



SH2208 First Aid Plaster

The product is low in sensitization, excellent in biological compatibility, bacteria-free and high in air permeation, and doesn't irritate or damage the skin nor affect the natural functions of the skin.
Material: PU, PVC, PU, Elastic Fabric, Spun-laced
Size: Ø22mm, 38×38mm, 40×10mm, 72×19mm, 76×19mm, 72×25mm, 76×30mm, 76×50mm, 76×38mm.

SH2301 Alcohol Prep Pad (Alcohol Swab)

Contents: 70% Isopropyl Alcohol or 75% Ethyl Alcohol
Non-Woven Size : 65×30mm, 65×58mm, etc
Type: 2ply, 4ply, etc.
Packet Size: 50×50mm
Packing: 100pcs/box, 200pcs/box, etc.



SH3102 Self-Sealing Flat Pouch

Size: 50×130mm, 90×230mm
135×255mm, 190×330mm
250×370mm



SH2651 Underpad

Material: non-woven fabrics+2ply medical paper + fluff pulp+PE film, non-woven fabrics+4ply medical paper+PE film, etc
with or without SAP
Size: 40×60cm, 60×90cm, etc
Weight can be customized.





SH2622 Disposable Isolation Gown

Material: SPP or PP+PE non-woven fabrics
or as customers' requirement
Color: Blue, yellow, green and so on.
Size: S/M/L/XL/XXL
Details: knitted or elastic cuff
With velcro or ties on the neck and ties on the waist.

SH2626 Disposable Protective Coverall

Chemical protective coverall includes: liquid-tight protective clothing , spray-tight protective clothing protective clothing for use against solid particulates protective clothing against liquid chemicals and so on.
Material: PP/PP+PE/SMS and so on
Weight: 35-60g
Size: M/L/XL/XXL
Color: white
Application: provide reliable protective products for industry/lab/ nuclear radiation/ agriculture/ military combat readiness and so on



SH2616 N95 Face Mask (with/without valve)

Material: polypropylene (conforms to American NEL SON N95/FFP1/FFP2/FFP3 test standards)
Specification: ear-loop
Size: without valve or with valve
Color: white



SH2615 Face Guard

Material: PET
Type: elastic
Double-sided anti-fog
Packing: 1pc/pe bag



SH0840 Virus Collection and Preservation System

Adapted to collection and transportation for Coronavirus, virus of flu, Bird flu, Hand-foot and mouth disease, measles etc.
Generally collect specimen from: Oral cavity, throat, nasopharynx, anus etc. System components:
1. Flocked swab: Oral type or Throat type; ABS Stick, Nylon Tip, 150mm
2. Virus preservation medium: Inactivated or Non-inactivated; 13x84mm or 16x100mm
3. Biohazard bag.



**SH0651
Tourniquet**

Highly flexible
Latex free cotton ribbon which is gentle on the skin
Sturdy snap-on clip
Quick and safe releasing



**SH06521
Non-latex Tourniquet**

Material: TPE
Size: 1"x18" , 0.6mm in thick, with/without texture
1pc/roll or Flat-Packing



**SH06522
Pre-cut Roll Tourniquet**

Material: TPE
Size: 1"x18" , 0.6mm in thick
25 each per big roll , per box



Pathology Grade
Non-interleaved
Dim.: 25.4mm x 76.2mm(1" x 3")
Thickness: 1.0 - 1.2mm

MICROSCOPE SLIDES
Color Frosted

Beveled Edges

Code	Description	Corner	
MF99138105W	White Color, Frosted one end one side	45°	
MF99138105O	Orange Color, Frosted one end one side	45°	
MF99138105Y	Yellow Color, Frosted one end one side	45°	
MF99138105G	Green Color, Frosted one end one side	45°	
MF99138105P	Pink Color, Frosted one end one side	45°	
MF99138105B	Blue Color, Frosted one end one side	45°	

Cut Edges

Code	Description	Corner	
MF99148105W	White Color, Frosted one end one side	90°	
MF99148105O	Orange Color, Frosted one end one side	90°	
MF99148105Y	Yellow Color, Frosted one end one side	90°	
MF99148105G	Green Color, Frosted one end one side	90°	
MF99148105P	Pink Color, Frosted one end one side	90°	
MF99148105B	Blue Color, Frosted one end one side	90°	



MF-TR-96 Filter Tips

Without any detectable DNA enzymes, RNA enzymes, proteins, DNA, heat and metal.

Size: 10uL, 20uL, 50uL, 100uL,

200uL, 1000uL, 1250uL

Aerosol-resistant

DNase/RNase Free

Filtered

Sterile

96 tips per rack

50 racks per case

- The tips are produced in Class 100,000 clean workshop, made from imported polypropylene material. Each batch is strictly evaluated, DNase-free, RNase-free, pyrogen-free, without endotoxin or PCR inhibitor.
- Universal Pyrogen-free pipette tip could be adapted to various pipettes like: Gilson, Eppendorf, Thermo (Finn), Keeven, Dragonlab, etc.
- The high quality pipette tip could simplify the strict workflow. Hydrophobic filter avoids cross-contamination of aerosols and liquids. Low-adhere tip could be adapted to viscose liquid and liquid with less surface tension, and reduce the liquid residue on the inner wall to ensure accurate research results.

Product Characteristics:




Better sealability with soft but tough centrum design.

Kinds of filters with different sizes of pore to ensure your result correct under different working situation.

Exact scale mark make reading easily and conveniently.

Tip go slender gradually ensure handling micro fluid.



Experiment prove tips with filter will work better under environment condition which with particle.



Vacuum Blood Collection Tube

CODE	TUBE NAME	CAP COLOR	TUBE SIZE(MM)	MATERIAL	ADDITIVE	DRAW VOLUME
LB6311B	PT TUBE (Double-deck)	BLUE	13x75	PP+PET	3.2% sodium citrate	1.8ML,2.7ML
LB6310R	PLAN TUBE	RED	13x75/13x100/16x100	GLASS/PET	NO ADDITIVE	1-10ML
LB6310O	PROCOAGULATION TUBE	RED/ORANGE	13x75/13x100/16x100	GLASS/PET	CLOT ADDITIVE	1-10ML
LB6310L	EDTA TUBE	LAVENDER	13x75/13x100/16x100	GLASS/PET	EDTAK2,EDTAK3	1-10ML
LB6310K	ESR TUBE	BLACK	13x75/8x120	GLASS/PET	3.8% sodium citrate	1.6ML,2.4ML 1.28ML,1.6ML
LB6310B	PT TUBE	BLUE	13x75/13x100	GLASS/PET	3.2% sodium citrate	1.8ML,2ML,2.7ML,3ML 3.6ML,4ML,4.5ML,5ML
LB6310H	GLUCOSE TUBE	GREY	13x75/13x100/16x100	GLASS/PET	POTASSIUM OXALATE +SODIUM FLUORIDE	1-10ML
LB6310G	HEPARIN TUBE	GREEN	13x75/13x100/16x100	GLASS/PET	SODIUM HEPARIN /LITHIUM HEPARIN	1-10ML
LB6310Y	GEL&CLOT ACTIVATOR TUBE	YELLOW	13x75/13x100/16x100	GLASS/PET	GEL+CLOT ACTICATOR	1-8.5ML

BLOOD COLLECTION NEEDLE



Multi-Sample Needle

The uniquely sharpened faceting simplifies penetration of tissue.
In comparison with other common needles, painful trauma is considerably reduced.
The specially sharpened point results in minimal tissue damage and thromboplastin seepage.

Code	Dsc.	Length	Color
MF6021 / MF6021S	21G	1"	Green
		1 1/4"	
		1 1/2"	
MF6022 / MF6022S	22G	1"	Black
		1 1/4"	
		1 1/2"	
MF6023 / MF6023S	23G	1"	Blue
		1 1/4"	
		1 1/2"	



Blood Collecttione Sets

Blood collection sets allows you to deliver exceptional patient care without compromising personal safety.
This product is intended for blood collection and infusion procedures. It is a necessity for automatic blood collection in the future.

Code	Dsc.	Length	Color
MF5921/MF5921S	21G	3/4"	Green
MF5922/MF5922S	22G		Black
MF5923/MF5923S	23G		Blue



MF5411

Mini Centrifuge

Max. Speed	7000rpm
Max. RCF	2680×g
Rotor Capacity	0.2/0.5/1.5/2.0mL×8; 0.2mL×16 PCR strips or 0.2mL×2 PCR 8 strips
Run Time	Continuous operation
Driving Motor	DC motor
Power	AC110-240V/50Hz/60Hz 20W
Noise Level	≤45dB
Dimension [D×W×H]	160x170x122mm
Weight	0.5kg



SH280109

Paper Safety Box

Capacity: 5L, 10L etc.



SH2802

Biohazard Specimen Bag

Size: 8"x10"
Thickness : 2mil 3ply
Size: 6"x10"
Thickness: 2mil 3ply